

## Product datasheet for **VC100067**

### E4 ORF6/7 (AC\_000007) Virus Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** E4 ORF6/7 (AC\_000007) Virus Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** E4 ORF6/7  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >The Viral ORF clone VC100067 represents NCBI reference of AP\_000191 with codon optimized for human cell expression  
Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCCGATCGCC

ATGACTACCAGCGGGGTGCCTTTCGGCATGACTCTCCGGCCTACCAGGTACGCTTGTACGCAGAACTC  
 CGTATTCACGGGACCGGCTCCCCCTTTGAGACCGAACTAGGGCCACAATTCTGGAGGATCACCCCTT  
 GCTGCCTGAATGTAACACGCTCAGATGCACAACGCCTGGACTTCTCCATCACCTCCAGTTGAACAACCA  
 CAGGTGGGCCAGCAGCCCGTGGCCCAACAACTTGATAGCGATATGAACTTGTCGGAAGTGCCTGGGAGT  
 TCATTAATATAACCGACGAAAGACTCGCCAGGCAGGAGACGGTCTGGAACATCACGCCAAAAACATGAG  
 CGTACTCACGATATGATGCTGTTTAAAGCTAGTCGCGGTGAGCGCACCGTATATTCTGTCTGCTGGGAG  
 GGTGGAGGGAGGCTCAACACACGGGTGCTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >VC100067 representing AP\_000191  
Red=Cloning sites Green=Tags

MTTSGVPGMTLRPTRSRLSRPTPYSRDLPPFETETRAILEDHPLLPECNTLTMHNAWTSPPVPEQP  
 QVGQQPVAQQLDSDMNLSELPGEFINITDERLARQETVWNIIPKNMSVTHDMMLFKASRGERTVYSVCWE  
 GGGRLNTRVL

TRTRRLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** AC\_000007

**ORF Size:** 450 bp

**OTI Disclaimer:** The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [AC\\_000007.1](#), [AP\\_000191](#)

**RefSeq ORF:** 450 bp

**UniProt ID:** [P03254](#)

**MW:** 17.1 kDa