

Product datasheet for VC100100

E3 RID-beta (AC_000019) Virus Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: E3 RID-beta (AC_000019) Virus Tagged ORF Clone
Tag: Myc-DDK
Symbol: E3 RID-beta
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >The Viral ORF clone VC100100 represents NCBI reference of AP_000598 with codon optimized for human cell expression
 Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCAGGCCATCTGCCATCTTCTCTTGCTGCTCCTCCCCTACGCCGATCTACCCCTGCAGCTTACT
 CTACTCCCCCGAACATCTTAGAAAATGCAAATCCAACAGCCCTGGAGCTTTCTGGCCTGCTACAGAGA
 AAAGTCGAAATCCCCCGAATTTGATCATGATCGCAGGCATCATTAAACATCATTGTGCTGTACCATTATC
 TCCTTTCTGATCATCCACTGTTTGACTTCGGTTGGAATGCTCCCAATGCCATGATCACCCACAAGATC
 CTGAAGAGCATATTCCACCGCAGAACATGCAGCACCCAATTGCTCTCATCGACTACGAAAGCGAGCCTCA
 GCCACCACTCCTGCCGCTATTTTCATACTTCAATCTTACGGGAGGGGATGAC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >VC100100 representing AP_000598
 Red=Cloning sites Green=Tags

MQAILPIFLLLLLPYAVSTPAAYSTPPEHLRKCKFQQPWSFLACYREKSEIPPNLIMIAGIINIICCTII
 SFLIYPLDFGWNAPNAHDHPQDPEEHIPPQNMQHPIALIDYSEPPPLLPASISYFNLTTGGDD

TRTRRLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: AC_000019

ORF Size: 402 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [AC_000019.1](#), [AP_000598](#)

RefSeq ORF: 402 bp

MW: 15.2 kDa