

Product datasheet for **VC100147**

Early E1A 12s (NC_012959) Virus Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Early E1A 12s (NC_012959) Virus Tagged ORF Clone
Tag: Myc-DDK
Symbol: Early E1A 12s
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >The Viral ORF clone VC100147 represents NCBI reference of YP_003038598 with codon optimized for human cell expression
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCGGCACCTCCACCTGCTGAGCTCTACTGTGCCTATCGATATGGCAGCCCTCCTTTGGAAGATTATG
TGAACACTATACTCGAGGACGAACCTCACTGTCTCCTTTTGAGCTCGGCCAACTCTGCAGGACCTGTA
TGATCTGGAAGTCAATGCCAGGAAAACGACCCCAATGAGGAGGCCGTTAACTTGATCTTCCCTGAGTCC
ATGATCTGCAGGCCGACATTGCAAGCGAAGCCGTTCTACATCTGTCTATACTCCCACTTTGCCACCCA
TTCCCAAGCTGGAGGAAGAGGATAAGCTTGATTTGAGGTGCTATGAGGAGGGATTCCACCGAGCGACAG
TGAAGATGAGAGAGGACCTGTGTCTGAAGACGAAAGTAGTCCCTGGAGGAGGACCATCCCTCCCTCCT
GATCTGACCAATGACACCCCTGCAAGTACGAAACCCACTCCTGTTAGACCAAGCGGAGAGCGACGCG
CCGCAAGTACAAAATTGAGGATCTGCTGCAGGACGTGGGGGAAACGAGCCCTGAACCTGTCTCTTAA
GAGACCCCGAAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >VC100147 representing YP_003038598
Red=Cloning sites **Green**=Tags

MRHLHLLSSTVPIDMAALLLEDYVNTILEDELHLSPFELGPTLQDLYDLEVNAQENDPNEEAVNLIFPES
MILQADIASEAVPTSVYTPPLPPIPKLEEDKLDLRCYEEGFPPSDSEDERGPVSEDESSPLEEDHPSPP
DLTNDTPLQVRKPTVPRPSGERRAAVYKIEDLLQDVGGNEPLNLSLKRPRN

TRTRRLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NC_012959

ORF Size: 573 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NC_012959.1](#), [YP_003038598](#)

RefSeq ORF: 573 bp

MW: 21.4 kDa