

Product datasheet for VC100160

pVI (NC_012959) Virus Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	pVI (NC_012959) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	pVI
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC100160 represents NCBI reference of YP_003038611 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGGAGGATATCAATTCGCCTCTCTGGCTCCACGGCATGGTACAAGGCCGTTTCATGGGGACTTGGAAATG
AGATCGGGACCTCTCAGCTTAACGGCGCGCTTTAACTGGAGCAGCGTATGGTCCGGACTTAAAAATTT
CGGGAGCACTTTGCGGACCTACGAAATAAAGCTTGAATAGCAGCACCAGGATGCTTAGAGAGAAA
CTCGTGACCAGAATTTCCAGCAGAAGGTCGTAGACGGGTTGGCCAGCGGATTAACGGTGTGTAGACA
TTGCCAACCGGCGTGCAACGAGAAATTAACCTACGCCCTGAATCCAAGGCCACCCACCGTGGTGGAGAT
GGAGGACGCCACACCCCTCCACGAGGGGAGAAACGCCCCCGGCCAGATGCCGAAGAAACGACCTTGCAG
GTAGACGAACCCCGAGCTACGAAGAGGCTCTCAAAGCTGGGATGCCTACCACGCGGATTATAGCCCAA
TGGCCACAGGCGTGATGAAACCCGCTACACTCGATCTCCATCTCCTCTGTTCTGCTCCTAAGGC
CGCCCCGGTGGTTCAGGCCACCCGTCGACAGCTGTGCGAAGGGTGCCCGCTAGGAGGCAAGCTCAG
AACTGGCAGTCCACCCTGACTCTATCGTGGGCCTTGGAGTTAAATCTCTGAAACGCCCGCGATGCTAT

ACGCGTACGGCGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >VC100160 representing YP_003038611
 Red=Cloning sites Green=Tags

MEDINFASLAPRHGTRPFMGTWNEIGTSQLNGGAFNWSSVWSGLKNFGSTLRTYGNKAWNSSTGQMLREK
 LADQNFQQKVVVDGLASGINGVVDIANQAVQREINSRLNRPPTVVEMEDATPPPGRKEKRRPDAEETTLQ
 VDEPPSYEEALKAGMPTTRIIAPMATGVMKPATLDLPSPVPAPPKAAPVVQAPPVATAVRRVPARRQAQ
 NWQSTLHSIVGLGVKSLKRRRCY

TRTRRLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NC_012959

ORF Size: 699 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NC_012959.1, YP_003038611
RefSeq ORF: 699 bp
MW: 25.4 kDa