

## Product datasheet for VC100190

### pX (AC\_000018) Virus Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** pX (AC\_000018) Virus Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** pX  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >The Viral ORF clone VC100190 represents NCBI reference of AP\_000546 with codon optimized for human cell expression  
 Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGCACTGACATGTAGACTGCGCGTGCCTATTACCGGATATCGAGGAAGAAATAGCAGGCGCCGGAGGA  
 TGCTGGGTGGGGTATGCGCCGCCACCGACGGAGACGCGCCATTTCCAGGAGGTTGGGGGCGGTTTTCT  
 CCCTGCTCTATCCCATCATCGCAGCTGCTATTGGGGCCATCCCTGGCATTGCTTCTGTGCGGGTGCAG  
 GCATCTCAACGCCAT

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >VC100190 representing AP\_000546  
 Red=Cloning sites Green=Tags

MALTCRLRVPITGYRGRNSRRRMLGRGMRRHRRRRAISRRLGGFLPALIPIIAAAIGAIPIGASVAVQ  
 ASQRH

**TR**TRRLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** Sgfl-MluI



[View online »](#)

Cloning Scheme:



ACCN: AC\_000018

ORF Size: 225 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [AC\\_000018.1](#), [AP\\_000546](#)

RefSeq ORF: 225 bp

MW: 8.3 kDa