

Product datasheet for **VC100219**

E1A (AC_000005) Virus Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	E1A (AC_000005) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	E1A
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC100219 represents NCBI reference of AP_000107 with codon optimized for human cell expression Red =Cloning site Blue =ORF Green =Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGAACCGAAATGACTCCCCTCGTTCTGTCCTACCAGGAGGCTGACGACATCCTGGAGCACTTGGTGG
ACAACCTCTTTAACGAGGTGCCTTCCGATGACGATCTGTACGTTCCCTCTCTATGAACTGTATGATCT
CGATGTGGAAAGCGCCGGCGAGGATAACAACGAACAGGCAGTAAATGAGTTTTTCCAGAGTCCCTTATC
TTGGCCGCGTCCGAAGGACTGTTTCTGCCGAGCCCCCGTCTTGAGCCCTGTGTGCGAACCCATTGGCG
GAGAATGTATGCCCCAGCTGCACCTGAGGATATGGATCTCTTGTGCTATGAAATGGGTTTTCCATGCTC
TGACAGCGAGGACGAACAGGATGAAAACGGAATGGCACATGTGTCCGCCAGTGCTGCAGCAGCTGCTGCT
GATCGGGAACGGGAGGAGTTCCAGCTGGACCACCCAGAGCTTCCAGGCCATAACTGCAAGAGCTGCGAGC
ACCACCGCAATCCACCGGTAACACAGATCTGATGTGTTCCCTTGTGTTATCTTAGAGCCTATAACATGTT
CATTTATAGCCCCGTTTCCAGACAACGAGCCAGAGCCCAATAGCACCTCGACGGAGATGAAAGGCCAAGC
CCTCCTAAGTTGGGATCCGCTGTACCCGAGGGAGTCATCAAACCCGTTCTCAGCGCGTGACCGGTGCA
GAAGGTGCGCCGTGGAGTCAATCTGGATTTGATTGAGGAAGAGGAAAGAGAACAGACAGTGCCGGTTGA
CCTTCCGTGAAAAGACCTCGATGCAAC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >VC100219 representing AP_000107
Red=Cloning sites Green=Tags

MRTEMTPLVLSYQEADDILEHLVDNFFNEVPSDDDLVYPSLYELYDLDVESAGEDNNEQAVNEFFPESLI
 LAASEGLFLPEPPVLSVPCPEIGGECMPQLHPEDMDLLCYEMGFPCSDSEDEQDENGMAHVSASAAAAA
 DREREEFQLDHPQLPHNCKSCEHHRNSTGNTDLMCSLCYL RAYNMF IYSPVSDNEPEPNSTLDGDERPS
 PPKLGS AVPEGVIKPVQRVTGRRRC AVESILDLIQEEREQTVPVDLSVKRPRCN

TRTRRLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: AC_000005

ORF Size: 798 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: AC_000005.1, AP_000107
RefSeq ORF: 798 bp
MW: 29.7 kDa