

## Product datasheet for VC100307

### L2 (NC\_011202) Virus Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	L2 (NC_011202) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	L2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC100307 represents NCBI reference of YP_002213808 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGTCCATTCTTATTTACCTTCTAACAAACAGGGTGGGGGCTGCGCGCCCATCAAAGATGTACGGAG  
GCGCACGAAAACGCTCTACCCAGCATCCCGTGCGGGTGCGGGGGCATTCCGCGCTCCCTGGGGTGCCCT  
CAAGGGACGGACACGAGTTCGCACAACCGTGGACGACGTGATCGATCAGGTCGTGGCAGACGCTCGGAAT  
TACACCCCAACCGCTCTACTTCAACAGTGGACGCGGTGATCGACTCAGTGGTGGCCGACGCACGCAAT  
ATGCAAGCGCAAGTCTCGGCGGAGACGAATTGCACGACGACACAGGGCTACAACAGCAATGAGGGCGGC  
CAGAGCTCTGCTGAGACGAGCTAGAAGAGTTGGGAGAAGAGCTATGCTGCGCGCTGCCAGACGCGCAGCA  
TCTGGTGCATCCGCCGGGCGGTCCCGCCGCCAGGCAGCAGCAGTTGCAGCTGCTACAATCGCTGATATGG  
CTCAAAGCCGACGAGGAAACGTTTATTGGGTCAGGGACGACGCAACCGCCAGCGAGTCCCGTCCGAAC  
TAGACCACCCCGCACT

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:	>VC100307 representing YP_002213808 Red=Cloning sites Green=Tags
-------------------	---

MSILISPSNNTGWGLRAPSKMYGGARKRSTQHPVVRVGRHFRAPWGALKGRTRVRTTVDDVIDQVVADARN  
YTPTAPTSTVDAVIDSVADARNYARRKSRRRRIARRHRATTAMRAARALLRRARRVGRRAMLRAARRAA  
SGASAGRSRRQAAAVAAATIADMAQSRRGNVYWRDAATGQRPVVRTRPPRT

**TR**TRRLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NC\_011202

**ORF Size:** 576 bp

**OTI Disclaimer:** The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NC\\_011202.1](#), [YP\\_002213808](#)

**RefSeq ORF:** 576 bp

**Locus ID:** 6870506

**MW:** 21.3 kDa