

## Product datasheet for VC100315

### E3 (NC\_011202) Virus Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** E3 (NC\_011202) Virus Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** E3  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >The Viral ORF clone VC100315 represents NCBI reference of YP\_002213825 with codon optimized for human cell expression  
 Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGCAAGCGATCCTTCCCATCTTTCTGCTGCTTTTCTGCTGCCATATGCCGTTTCCACACCTGCAGCGTATT  
 CAACCCCTCTGAACATCTTAGAAAATGCAAGTTTCAACAGCCCTGGTCTTTTCTGGCGTGCTATCGCGA  
 AAAATCCGAAATCCACCAAACCTGATTATGATAGCAGGGATTATCAACATTATTTGCTGCACAATTATA  
 AGCTTTCTGATTACCCCTGTTTCGATTTTGGATGGAATGCACCGAATGCCACGACCACCTCAGGATC  
 CAGAGGAACACATTCCGCTCCAGAATATGCAACACCCCATCGCCCTCATCGATTATGAATCCGAGCCACA  
 GCCTCCTCTGCTGCCGCCATTAGCTACTTCAATTTGACCGTGGAGATGAC

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >VC100315 representing YP\_002213825  
 Red=Cloning sites Green=Tags


MQAILPIFLLLLLPYAVSTPAAYSTPPEHLRKCKFQQPWSFLACYREKSEIPPNIIMIAGIINIICCTII  
 SFLIYPLDFGWNAPNAHDHPQDPEEHIPLQNMQHPIALIDYESEPPPLPAISYFNLGGDD

**TR**TRRLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** Sgfl-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NC\_011202

**ORF Size:** 402 bp

**OTI Disclaimer:** The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NC\\_011202.1](#), [YP\\_002213825](#)

**RefSeq ORF:** 402 bp

**Locus ID:** 6870511

**MW:** 15.2 kDa