

## Product datasheet for **VC100351**

### E3 (NC\_011203) Virus Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** E3 (NC\_011203) Virus Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** E3  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >The Viral ORF clone VC100351 represents NCBI reference of YP\_002213787 with codon optimized for human cell expression  
 Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGAAAGCTTTCGCCGTGCTGTTCTGACTGAGCCTCATCAAACGGAAGTGGACCCCTCATATGGGCTGC  
CCCTCTCAAAGCGGTTTGTACAATACGAATCAGATTTCCAGAAAACCCAGACCCTGCCCCCTTCAT  
TCAAGACTCAAATTCTACACTGCCCGCACCCCTCAACAACCAATCTGCCAGAAACCAATAAGTTGGCAAGT  
CACCTGCAGCACAGGCTCTCCAGGAGTCTGCTGAGTGCCAATACCACTACTCCTCAAAAACAGGCGCGCAAC  
TGAGAGGTTTGCCGACGGATGACCCATGGGTTGTGCGCGGATTTCGTAACGCTGGGCGTCGTGGCTGGAGG  
CTTGGTGCTGATACTTTGCTACCTGTACTCCGTGTTGCGCATATCTGGTGATCCTCTGCTGCTGGTTC  
AAAAAGTGGGACCATAT

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >VC100351 representing YP\_002213787  
 Red=Cloning sites Green=Tags

MKAFAVLFVLSLIKTELPSYGLPLLQSGLYNTNQIFQKQTQLPPFIQDSNSTLPAPSTTNLPETNKLAS  
HLQHRLSRSLLSANTTTPKTGGELRGLPTDDPWVAVGFVTLGVVAGGLVILCYLYTPCCAYLVILCCWF  
KKWGPY

**TR**TRRLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NC\_011203

**ORF Size:** 438 bp

**OTI Disclaimer:** The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NC\\_011203.1](#), [YP\\_002213787](#)

**RefSeq ORF:** 438 bp

**Locus ID:** 6870528

**MW:** 16.0 kDa