

Product datasheet for **VC100443**

L4 (NC_010956) Virus Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: L4 (NC_010956) Virus Tagged ORF Clone
Tag: Myc-DDK
Symbol: L4
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >The Viral ORF clone VC100443 represents NCBI reference of YP_001974455 with codon optimized for human cell expression
Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCCGATCGCC

ATGCCAAGGAAGAAGCAAGAGCCTCTGGTGAAGAAATGGAAGAAGAATGGGACAGTCAAGCTGAGGAGGATGAATGGGAGGAAGAACGGAGGAAGAAGAGTTGGAAGAGGTGGAGGAGGAGCAAGCTACAGAACAGCCCGTCGCAGCTCCATCTGCACCCGCCGCCACAGCCGTCACCGATACAACACTAGTGCTGCGCCGCGAAGCCCAAGGAGGTGGGATCGGGTGAAAGGAGATGCTAAAAAAGCAGGTGCGGGAGGTGCGCCGGGGTGGCCTGAGGATCGCCGCAAACGAACCCCTCACAAACACGCGAGCTTCGCAATCGAATTTTCCGACCCTTTATGCATCTTCCAGCAGTCCCAGGCCAACAGCAAGAGCTGAAAGTTAAGAATCGCAGTCTCCGGTCCCTCACAGGAGCTGTGTACCATAAGAATGAGGACCAATTGCAGCGCACCTCGAGGATGCCGAGCGCTGTTCATAAATATTGCGCTTTGACTCTGAAGGAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >VC100443 representing YP_001974455
Red=Cloning sites Green=Tags

MPRKKQEPLVEEMEEWDSQAEDEWEETEEEELEEVVEEQATEQPVAAPSAPAAVTDTTSAAPAKP
 PRRWDRVKGDAKKKQVREVAGGGLRIAANEPLTTRELRNRIIPTLYAIFQQSRGQQQELKVNRSLSL
 RSCLYHKNEQLQRTLEDAEALFHKYCALTLKD

TRTRRLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:


ACCN: NC_010956

ORF Size: 519 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NC_010956.1](#), [YP_001974455](#)

RefSeq ORF: 519 bp

Locus ID: 6386277

MW: 19.9 kDa