

Product datasheet for **VC100447**

E3 (NC_010956) Virus Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	E3 (NC_010956) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	E3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC100447 represents NCBI reference of YP_001974458 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCGCATTTTTGCTGTCCTGTTCTGTGGTATCACTGATTAAGCTAAACTGTTGCAATACTCTGGCCCAT
GCCGGCACCAGTTGAATGAGACCGTTTACCTGACCAACCAAACCTGAGGTTAAGCTCACTTGCCGCCACACA
CAAGACCTACATAATCTGGTTTTTCGAGAACACTTCATTTGCCGTAGCAACACACATTGCAATGACGGC
GTGGAGCTCCGAACAATCTGAGCTCAGGCCTTCCTATGACACTCATAGAGCGAAACTCGTGCTGTATA
ACCCCTTTGTCGAGGGGACCTATCACTGTCACTGTGAGTCTGGGCCGTGCTTCCATTCTTTACACTGGTGAACGT
GACAGGATCTTCAACCGCCGCCCCGAGCCCAATCTCCCTAGTGACACAATCAAACCAAGATTTGGCGGA
GAGCTGAGACTGCCACCTAGCGAAGAAGGGGTTTCTCCATATGAAGTGGTGGGTACCTGATCCTCGGAG
TCGTCCTTGGTGGCTGTATCGCGGTCCTGGCCAGTTGCCCTGTTGGGTTGAGATTAATAATTTTTATTTG
TTGGGTGAGACTGCGGAGAGGAGCTG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:	>VC100447 representing YP_001974458 Red=Cloning sites Green=Tags
-------------------	---

MRIFAVL FVVS L I KAKLLQYSGPCRHLNETVYLTNQTEVKLTCRPHKTYIIWFFENTSFVSNTHCNDG
VELPNNLSSGLSYDTHRAKLVLYNPFVEGTYHCQSGPCFHSFTLVNVTGSSTAAPENLPSDTIKPRFGG
ELRLPSEEGVSPYEVVGYLILGVVLGGCIAVLAQLPCWVEIKIFICWVRHCGEEL


TRTRRLEQKLI SEEDLAANDILDYKDDDDKV



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NC_010956

ORF Size: 588 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NC_010956.1](#), [YP_001974458](#)

RefSeq ORF: 588 bp

Locus ID: 6386273

MW: 21.8 kDa