

Product datasheet for **VC100626**

3C (NC_012801) Virus Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	3C (NC_012801) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	3C
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC100626 represents NCBI reference of YP_002956114 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGCCCGATCTGCGAGACGTGGTGACTTTCTTCAACGCAGCCAAAGGAGCCAGTGGATGGTGGAGT
CCATCAAGAGCCTGATAAACTGGATTAAGCAGTGGCTGGAGTTGGAAGAGGAGAACGAGGCAGTTAAGCT
GGAGAAAATGCTGATCGAGTCACCATCACACTGTAAAAACATACACATGTACAACAAGGGCGAATTGTTC
GCCAAACCTCAAGAGTCTTTGACTTCTGGACCGGCTGTGCGAGACAGCTACAAGCCTCGGTAAGACAC
ACCTGGCATCTTACTTCCGCAATTTCTGTGACATACGACTCAGATACATCTCGCCCTGAACCAAGTCGTAGT
TGTGTTGAGGGTAAGCCCGGAGCCGGCAAATCAGCCGCTGCCACCGTCTTGGCTGCAGCAGTATCAAAG
CTGTGTGTCGGCAGTCAATCAGTCTACACATTGTCCCCGACACAGAACATATGGACGGCTATCACGGGC
AGTTCGCAACCATCATGGATGATCTGGCCAGAACCCTGACGGTGAAGATTTCCGCTCTTTCTGTGAGT
GGTTAGCTGCGCCAGTACCGGCCGCTATGGCAGGTCTCCAGATAAGGGTATACTCTTACCTCCCGG
ATTGTGATCGCGACCACCAATCTTTGCGACTTTAACCCTTGACTATAGCCGACCCTAGAGCCCTGGAGC
GGCGCATCACCTTTGATATTCTGGTCTCACCAGGTCTGCCTGCAACAAAATGGAAAACCGACCTGAA
CGCTGCCCTGACACCTGATGGACCAGGTGAAGGGCCATTCAACTGATTGCCAAATTCGCACACAACA
GGACTGACTCTGAAGGATATCCGCAAAGGAACCACTGAACCTCAAGGACCTGACAGAAGTGGTGGTTG
ATAGAATTAACAAGAAGAAGAAGGTGGGAACATGCTTGAGAACCTGGTCTGCTCAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >VC100626 representing YP_002956114
Red=Cloning sites Green=Tags

MGPLDRDVTFFNAAKGAQWMVESIKSLINWIKQWLEEEENEAVKLEKMLIESPSHCKNIHMYNKGELF
 AKPQESFDLDRLCETATSLGKTHLASYFRNFVTYSDTSRPEPVVVVLRGKPGAGKSAATVLA AVSK
 LCVGSQSVYTLSPDTEHMDGYHGQFATIMDDLQNPDGEDFRSFCQMVSCAQYRPPMAGLPDKGILFTSR
 IVIATTNLCDFNPLTIADPRALERRITFDILVSPGPACNKNGLDLNAALTPDGPGEPPFTTDCQILHTT
 GLTLKDIRKGTTLNLKDLTELVDRINKKKKVGNMLENLVAQ

TRTRRLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NC_012801

ORF Size: 966 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NC_012801.1, YP_002956114
RefSeq ORF: 966 bp
MW: 35.5 kDa