

Product datasheet for **VC100661**

HECV-4408 4.8 kD NSP Gene Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	HECV-4408 4.8 kD NSP Gene Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	4.8 kD non-structural protein
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC100661 represents NCBI reference of YP_003038524 with codon optimized for human cell expression Red =Cloning site Blue =ORF Green =Tags(s) GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCC CGATCGCC ATGCCAATGGCAACA ACTATAGATGGAACCGATTATACTAATATCATGCCCTCTACTGTGTTCACTACGGTGTATCTCGGTGGCTTTATTGGCATTGACACTTCAACGACCGGCTTTACTTGTTTTCAAGGTAC ACGCGT ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
Protein Sequence:	>VC100661 representing YP_003038524 Red =Cloning sites Green =Tags MPMATTIDGTDYTNIMPSTVFTTVYLLGGFIGIDTSTTGFTCFSTRY TRTRRLEQKLI SEEDLAANDILDYKDDDDKV
Restriction Sites:	Sgfl-MluI



[View online »](#)

Cloning Scheme:



ACCN: NC_012950

ORF Size: 135 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NC_012950.1](#), [YP_003038524](#)

RefSeq ORF: 135 bp

MW: 4.9 kDa