

## Product datasheet for **VC100713**

### 3C (NC\_013114) Virus Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** 3C (NC\_013114) Virus Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** 3C  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >The Viral ORF clone VC100713 represents NCBI reference of YP\_003104784 with codon optimized for human cell expression  
 Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGGACCTGCTTTCGAGTTCGCCGTAGCCATGATGAAACGCAACGCAAGTACTGTGAAGACGGAGTACG  
GTGAATCACCATGCTTGGGATTTACGATAGGTGGCAGTACTCCAAGGCATGCCAAACCAGGACCCAC  
GATCCTGATGAATGATCAGGAGGTCGGTGTCTGGATGCAAAGGAGCTCGTGGATAAGGACGGGACAAC  
CTGGAACCTTACACTGCTGAACTGAACAGGAATGAAAAGTTTAGGACATCAGGGGATTCCTGGCCAGAG  
AGGAAGCCGAAGTGAATGAGGCCGTAATCGCAATCAATACTTCCAAGTCCCAACATGTATATTCTGT  
GGGGCAAGTCGCCGATTACGGTTTCCTGAATCTGGGCGGACCCCAACCAAGCGGATGCTGATGTACAAC  
TTTCCAACACGCGCCGTCAGTGTGGTGGCGTCTCATGAGTACAGGAAAGGTCCTGGGCATCCACGTTG  
GCGGCAATGGACATCAGGGATTCTCCGAGCTTTGCTTAGGCATTACTTCAACGATGAGCAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >VC100713 representing YP\_003104784  
 Red=Cloning sites Green=Tags

MGPAFEFVAMMKRNASTVKTEYGEFTMLGIYDRWAVLPRHAKPGPTILMNDQEVGVVDAKELVDKDGTN  
LEL TLLKLNREKFRDIRGFLAREEAEVNEAVLAINTSKFPNMYIPVQVADYGFNLGGTPTKRMLMYN  
FPTRAGQCGVLMSTGKVLGIHVGGNGHQGFSALLRHVFNDEQ

**TR**TRRLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI



**Cloning Scheme:**


**ACCN:** NC\_013114

**ORF Size:** 552 bp

**OTI Disclaimer:** The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NC\\_013114.1](#), [YP\\_003104784](#)

**RefSeq ORF:** 552 bp

**MW:** 20.4 kDa