

## Product datasheet for **VC100767**

### capsid protein 1C version 1 (NC\_001430) Virus Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	capsid protein 1C version 1 (NC_001430) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	capsid protein 1C version 1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC100767 represents NCBI reference of NP_740739 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGGGTTCTACTTGTCTGCTGCCCGGAGCAATCAGTTCCTGACCACTGACGACCACAGTTCCTGCC  
CAGCTTTCCAGATTTTAGTCCCACACCTGAGATGCACATCCCTGGACAGGTCCATTCCATGCTCGAGAT  
AGTCCAGATTGAGAGCATGATGGAAATTAATAACGTGAATGACGCCTCAGGGGTGGAACGATTGCGGGTC  
CAGATATCAGCCAGTCAGACATGGATCAGCTCCTGTTAATATCCCTCTTGATATACAGCTCGAGGGTC  
CTTTGCGAAATACACTGCTGGGCAACATTAGCCGCTACTACACTCACTGGAGCGGGTCCCTGGAAATGAC  
ATTCATGTTCTGCGCAGCTTCATGACTACCGGTAACACTGATAATTTGCTATACCCCTCCCGGGGCTCA  
TCCCCAACTGATAGGATGCAGGCAATGCTGGCTACACACGTGGTCTGGGATTTTGGGCTTCAGTCTCAA  
TCACCATTATAATCCCCTGGATCAGTGGGAGCCATTACCGCATGTTAACACCGATGCTAAGGCTATTAA  
TGCTAACGTAGGCTATGTTACCTGCTTCATGCAAATAACCTCGTCGCGCCCGTTGGCGCTGCAGACCAA  
TGTTATATTGTTGGAATGGTTGCAGCTAAGAAGGACTTTAATCTGAGACTTATGCGGGATTCCCCTGATA  
TCGGGCAGAGCGCAATTCTGCCGAACAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >VC100767 representing NP\_740739  
 Red=Cloning sites Green=Tags

MGVPTCLLPGSNQFLTDDHSSAPAFPDFSPTPEMHIPGQVHSMLEIVQIESMMEINNVNDASGVERLRV  
 QISAQSDMDQLLFNIPLDIQLEGPLRNTLLGNISRYYTHWSGSLEMTFMFCGSFMTTGKLIICYTPPGGS  
 SPTDRMQAMLATHVWDFGLQSSITIIIPWISGSHYRMFNTAKAINANVGYYTCFMQTNLVAPVGAADQ  
 CYIVGMVAAKKDFNLRMRDSDIGQSAILPEQ

TRTRRLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NC\_001430

**ORF Size:** 729 bp

**OTI Disclaimer:** The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NC\\_001430.1](#), [NP\\_740739](#)  
RefSeq ORF: 729 bp  
MW: 26.8 kDa