

## Product datasheet for **VC100775**

### protein 3A (NC\_001430) Virus Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** protein 3A (NC\_001430) Virus Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** protein 3A  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >The Viral ORF clone VC100775 represents NCBI reference of NP\_740419 with codon optimized for human cell expression  
**Red**=Cloning site **Blue**=ORF **Green**=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGC**C

ATGGGGCCACCAACTTTTAAGGAGATCAAATATCCGTCACCCCCGAACTCCTGCTCCTGATGCAATCA  
ACGATTTGTTGCGGTCCATCGACTCCCAGGAGGTACGCGATTACTGCCAAAAAAGGGGTGGATTGTGAT  
GCACCCGCCACGGAATTGTCGTCGATAAGCACATCTCAAGAGCCTTCATTGCATTGCAGGCTATAACG  
ACATTTGTGAGCATTGCCGGCGTGGTCTATGTGATCTACAAGCTGTTGCTGGTATACAA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >VC100775 representing NP\_740419  
**Red**=Cloning sites **Green**=Tags

MGPPTFKEIKISVTPETPAPDAINDLLRSIDSQEVRYDQCQKKGWIVMHPPTLVVDKHISRAFIALQAIT  
TFVSIAGVVVYIYKLFAGIQ


**TR**TRRLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI



[View online »](#)

Cloning Scheme:



ACCN: NC\_001430

ORF Size: 270 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NC\\_001430.1](#), [NP\\_740419](#)

RefSeq ORF: 270 bp

MW: 10.0 kDa