

Product datasheet for **VC100819**

UL31 (NC_001806) Virus Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	UL31 (NC_001806) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	UL31
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC100819 represents NCBI reference of NP_044633 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTACGATACCGACCCACATAGAAGGGGGTCTCGCCCCGACCGTATCATGGTAAGGAACGGCGGAGAT
CAAGGTCAAGCGCTGCAGGTGGTACTTTGGGAGTCGTGAGAAGAGCCAGTCGGAAAAGTCTTCCACCCCA
CGCTAGAAAACAGGAGCTCTGTCTGCATGAGCGCCAGCGCTATCGCGGCCTTTTGGCCGACTGGCTCAG
ACCCCAAGTGAGGAAATCGAATCGTCAGATCCCTTAGCGTGCCACTCGTCAAGACTACCCCTGTCAGCC
TGCCGTTCTGTCTCGATCAAACGGTGGCCGATAATTGCCTGACCCCTCAGCGGAATGGGGTATTACCTGGG
GATTGGAGGGTGCTGCCCTGCTTGAATGCCGGAGATGGCGGTTTGCCGCCACTTCAAGAGAGGCCCTC
ATCCTCGCTTCGTACAACAAATTAACACTATTTTGAACACAGAGCCTTTCTTGCGAGTCTTGTGGTCC
TTGCCGACCGGCACAACGCGCCACTTCAGGATCTCCTCGCCGGCATCCTCGGCCAGCCCGAACTGTTTTT
TGTACACACCATACTGCGAGGAGGCGGGCTTGCATCCTCGGCTGTGTTCTACCCTGACCCAACATAT
GGGGGCCACATGCTCTACGTAATCTTCCCGGCACTTCTGCCACTTGCATTATAGACTCATCGACAGGA
TGTTGACCGCTTGCCAGGATACAGATTCGTGGCCATGTCTGGCAGTCCACCTTCGACTGGTAGTGCC
ACGAAATGCCGAAAACCTACCGATGCAGAAATCCCCACCGTGTCCGCTGCAGACATTTACTGTAATG
AGGGACATCTCTTTGACGGCGGACTGATGCTCGAGTACCAAAGGCTGTATGCCACCTTCGACGAATTC
CACCTCCC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >VC100819 representing NP_044633
 Red=Cloning sites Green=Tags


MYDTPHRRGSRPGPYHGKERRRSRSSAAGGTLGVVRRASRKSLPPHARKQELCLHERQRYRGLFAALAQ
 TPSEEIAIVRSLSVPLVKTTVPVSLPFCLDQTVADNCLTLSGMGYYLGIGGCCPACNAGDRFAATSREAL
 ILAFVQQINTIFEHRAFLASLVVLADRHNAFLQDLLAGILGQPELFFVHTILRGGGACDPRLLFYPDPTY
 GGHMLYVIFPGTSAHLHYRLIDRMLTACPGYRFVAHVWQSTFVLVRRNAEKPTDAEIPVSAADIYCKM
 RDISFDGGLMLEYQRLYATFDEFPPP

TRTRRLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NC_001806

ORF Size: 918 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NC_001806.1](#), [NP_044633](#)

RefSeq ORF: 918 bp

Locus ID: 2703350

UniProt ID: [P04413](#)

MW: 34.0 kDa