

## Product datasheet for **VC100844**

### UL55 (NC\_001806) Virus Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	UL55 (NC_001806) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	UL55
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC100844 represents NCBI reference of NP_044658 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGACCGCCACGCCACTGACCAATCTTTTCTGCGCGCACCTGACATTACTCATGTGCCCCCCTTATT  
GCTTGAATGCTACTTGGCAAGCCGAAACCGCAATGCACACTTCTAAAACCGACTCCGCCTGTGTGGCCGT  
TAGAAGCTACCTGGTTAGGGCTAGCTGCGAGACAAGTGAACAATTCATTGTTTTTTTTGCGGTAT  
AAGGATACACATCACACTCCGCCTGATTACCGAATTGCGAACTTCGCCGACCTGTAAACCCAC  
CTGTGTTGAGGAATTGGAAGATAAGAGGGGGTAAGACTGAGGTGTGCGCGGCCCTTTTCCGTTGGGAC  
TATAAAGGACGTGTCGGCAGCGGAGCTTCTTACGCCGGGAGTATACGATTAATGGAATCGTCTATCAT  
TGCCATTGTAGATACCCATTCTCTAAAACCTGTTGGATGGGCGCCAGCGCCGACTGCAGCATCTCCGAA  
GTATCAGCAGCTCCGGCATGGCTGCTAGAGCTGCCGAGCATCGGAGGGTAAAAATTTAAATCAAGGCC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:	>VC100844 representing NP_044658 Red=Cloning sites Green=Tags
-------------------	--

MTATPLTNLFLRAPDITHVAPPYCLNATWQAETAMHTSKTDSACVAVRSYLVRASCETSGTIHCCFFAVY  
KDTHTTPLLITELRNFADLVNHPVLRLELDKRGVRLRCARPFVSGTIKDVSGSGASSAGEYTINGIVYH  
CHCRYPFSKTCWMGASAALQHLRSISSGMAARAAEHRRVKIKIKA

**TR**TRRLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:	SgfI-MluI
--------------------	-----------



[View online >](#)

**Cloning Scheme:**


ACCN: NC\_001806

ORF Size: 558 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NC\\_001806.1](#), [NP\\_044658](#)

RefSeq ORF: 558 bp

Locus ID: 2703427

UniProt ID: [P04413](#)

MW: 20.5 kDa