

Product datasheet for **VC100896**

UL31 (NC_001798) Virus Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	UL31 (NC_001798) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	UL31
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC100896 represents NCBI reference of NP_044501 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTATGACATTGCCCCCGCGCTCTGGATCTCGCCCTGGCCCCGAAGAGACAAGACTCGGAGGAGAT
CTAGGTTTTAGCCGCAGGCAACCCTGGCGTGGAAAGACGGGCCCTCCAGGAAGAGCCTGCCGAGTCATGC
CCGGCGACTGAACTGTGCCTCCACGAGCGCGACGCTATAGAGGCTTCTCGCTGCACTCGCACAGACA
CCATCCGAGGAAATTGCCATCGTGCGGAGCCTGTCTGTGCCTCTGGTAAAACACACCTGTCAGCCTCC
CTTTTAGTCTCGATCAGACGGTGGCAGATAAATTGTCTGACTCTCAGTGGTATGGGATACTATTTGGGCAT
CGGAGGCTGCTGCCCGCGTGTAGCGCCGAGACGGACGCTGGCTACTGTGTCTCGCGAGGCACTGATC
CTGGCCTTTGTCCAGCAGATAAAACCATTTCGAGCATAGGACTTTTCTCGCATCATTGGTTGTCTTGG
CTGACAGACACAGCACACCACTGCAGGACTTGCTTGTGATACCCTGGGCAACCCGAACTGTTTTTGT
GCATACTATCCTGCGGGAGGAGGCGCCTGCGACCCAGATTCTGTTTTATCCCGACCTACCTACGGG
GGACACATGCTGTACGTAATATTTCTGGGACCAGTGCGCATCTTCACTACAGACTTATTGACCGCATGC
TCACAGCATGTCCAGTTACCGTTTGCCGCACATGTCTGGCAGTCTACGTTTGTCTCGTCGTGAGGCC
GAATGCAGAGAAACCCGCGATGCCGAGATCCCAACAGTTAGCGCAGCCGACATATACTGTAATAATGCGG
GATATTAGCTTTGATGGCGGGCTGATGCTCGAATACCAGCGACTCTATGCTACATTGACGAATCCCC
CTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >VC100896 representing NP_044501
Red=Cloning sites Green=Tags

MYDIAPRRSGSRPGPRDKTRRRSRFSAAGNPGVERRASRKSLSHARRLELCLHERRRYRGFFAALAQT
 PSEEIAIVRSLVPLVKTPVSLPFLSLDQTVADNCLTLSGMGYYLGIGGCCPACSAAGDGRLATVSREALI
 LAFVQQINTIFEHRTFLASLVLADRHSTPLQDLLADTLGQPELFFVHTILRGGGACDPRFLFYDPPTYG
 GHMLYVIFPGTSAHLHYRLIDRMLTACPGYRFAAHVWQSTFVLVVRRAEKPADAEIPTVSAADIYCKMR
 DISFDGGLMLEYQRLYATFDEFPPP

TRTRRLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NC_001798

ORF Size: 915 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NC_001798.1, NP_044501
RefSeq ORF: 915 bp
Locus ID: 1487317
MW: 33.9 kDa