

Product datasheet for **VC100968**

ORF27 (NC_001348) Virus Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ORF27 (NC_001348) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ORF27
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC100968 represents NCBI reference of NP_040150 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGGATCTCTGCTTCGACGACATGCAGTACCGATGGAGCCCCTCCAATACACCCTACCGCTCAATGA
GTCGCCGCTACAAAAGTGTGCCCGAGCGGCCCTAGTATGAGAGTTCGGTCCAGGACTCCCTGCAGACG
CCAGACTATCCCGGAAAAGTGTGAGCAAAGAACGCAGTGTCTACAGACATTACTTTAATTACATCGCT
AGAAGTCCACCCGAGGAGCTTGCTACCGTCCGAGGCCGATTGTGCCGATTATAAAAACGACTCCAGTTA
CTCTGCCATTTAACTTGGGCCAGACTGTCGCCGACAACGTCTCTCACTCTCCGGGATGGGATATCATCT
TGGTCTCGGAGGGTATTGTCCAACCTGTACCGCCTCTGGCGAACCCCGCTTGTGTAGGACTGATAGAGCC
GCTCTGATCCTGGCCTATGTCCAGCAGCTCAATAATATATACGAGTACAGAGTATTTTTGGCCTCTATTC
TGGCGTTGTCCGACAGAGCTAACATGCAGGCTGCCTCCGCTGAGCCCCTTCTTTCTAGTGTGTTGGCCCA
GCCTGAACGTCTTTCATGTACCATATCATGCGGGAGGGCGGGATGCGCGATATCAGGGTCTTTTCTAC
AGGGACGGCGACGCGGGAGGCTTCATGATGTATGTCATTTTTCCAGGGAAGAGCGTTACCTCCATTATC
GCCTCATTGACCACATTCAGGCTGCGTGTGAGGATATAAGATTGTGGCTCACGTGTGCCAGACAACCTT
TCTTCTGTCAGTTGTGCAACCCTGAACAACAGACAGAGACCGTGGTCCCAGCATTGGGACAAGCGAT
GTCTACTGTAAGATGTGTGATCTTAACTTTGACGGTGAGCTGTTGCTCGAATACAAGAGACTGTATGCAC
TTTTCGACGATTTTGTGCCCCCAGG

ACGCGTACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >VC100968 representing NP_040150
Red=Cloning sites Green=Tags

MEDLCFDDMQYRWSPNTPYRSMSRRYKSVSRSGPSMRVRSRTPCRRQTIRGKLSKERSVYRHYFNIA
 RSPPEELATVRGLIVPIIKTTPVTLPFNLGQTVADNCLSLSGMGYHLGLGGYCPTCTASGEPRLCRTDRA
 ALILAYVQQLNNIYEYRVLASILALSDRANMQAASAEPLLSVLAQPELFFMYHIMREGGMRDIRVLFY
 RDGDAGGFMMYVIFPGKSVHLHYRLIDHIQAACRGYKIVAHVWQTTFLLSVCRNPEQQTETVVPSIGTSD
 VYCKMCDLNFDELLLEKRYLALFDDFVPPR

TRTRRLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NC_001348

ORF Size: 936 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NC_001348.1](#), [NP_040150](#)
RefSeq ORF: 936 bp
Locus ID: 1487666
UniProt ID: [P09310](#)

MW: 35.7 kDa