

Product datasheet for **VC101132**

BZLF2 (NC_007605) Virus Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	BZLF2 (NC_007605) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	BZLF2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC101132 represents NCBI reference of YP_401672 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGTGTCTTTAAGCAGGTTAGGGTGCCACTGTTTACCGCCATAGCACTGGTTATTGTGCTCCTGCTCG
CCTACTTCCTGCCTCAAGGGTGCAGGGAGGTGGCCGCGTGGCAGCTGCCGCAATTACATGGGTTCCCAA
ACCAAACGTGGAGGTGTGGCCCGTCGATCCACCACCGCCTGTGAATTTCAACAAAACGGCCGAGCAGGAA
TACGGAGACAAGGAGGTGAAGCTCCCCATTGGACTCCTACCCTTCATACATTTTCAGGTGCCCCAGA
ACTACACAAAAGCTAACTGCACGTATTGCAACACTAGGGAGTATACATTTAGTTACAAAGGATGTTGTTTTTA
TTTACCACAAAAAACAACACATGGAATGGATGCTTTCAAGCATGCGCTGAGCTGTACCCCTGCACCTAT
TTTTACGGTCTACCCCTGATATCCTGCCGTTGTGACACGAACTTGAACGCAATTGAGTCTCTGTGGG
TTGGGGTGTACCGGGTGGCGAGGGAACTGGACATCACTGGACGGGGAAACATTTAAAGTTTATCAGAT
ATTTGGCTCCCATTCACATACGTCTCCAAGTTCTCTACGGTCCCAGTGAGTCATCATGAATGTTCTTT
TTGAAACCCTGCCTCTGCGTTTCCAGAGGTCTAACAGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >VC101132 representing YP_401672
 Red=Cloning sites Green=Tags

MVSFKQVRVPLFTAIALVIVLLLAYFLPPRVGGGRVAAAAITWVPKPNVEVWPVDDPPPPVNFNKTAEQE
 YGDKEVKLPHWPTLHTFQVPQNYTKANCTYCNTRYTF SYKGCCFYFTKKKHTWNGCFQACAELYPCTY
 FYGPTDILPVVTRNLNAIESLVWGVYRVGEGNWTSLDGGTFKVVYQIFGSHCTYVSKFSTVPVSHHECSF
 LKPCLCVSQRSNS

TRTRRLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NC_007605

ORF Size: 669 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NC_007605.1](#), [YP_401672](#)
RefSeq ORF: 669 bp
Locus ID: 3783745
UniProt ID: [P03211](#)
MW: 25.3 kDa