

Product datasheet for **VC101170**

RPMS1 (NC_007605) Virus Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: RPMS1 (NC_007605) Virus Tagged ORF Clone
Tag: Myc-DDK
Symbol: RPMS1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >The Viral ORF clone VC101170 represents NCBI reference of YP_401710 with codon optimized for human cell expression
 Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCCGGGGCACGGAGACGCGCTAGATGCCAGCTAGTGCAGGATGTGCGTATAGCGCTCGACCTCCTC
 CGCTTAGCACCCGGGGCAGGAGAATCTCAGCTGGATCTGGCCAGCCAGGTGGTGGCCTTGGGGATCTCC
 TCCCCACCAGATACCCGCTACCGGCGACCTGGTCCCAGGAGAGAGCTAGGTCATGTCTGCACGCAGGG
 CCTCGCGGGAGACCGCCACACTCCCGAACAAGAGCTAGGAGGACCTCACCCGGGGCCGGAGGAGGAGGAT
 GGAGAGCGGGGAGCTGCACTTCCAGAGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >VC101170 representing YP_401710
 Red=Cloning sites Green=Tags

MAGARRRARCPCASAGCAYSARPPPLSTRGRRISAGSGQPRWWPWGSPPPDTRYRRPGRRRARSCLHAG
 PRGRPPHSRTRARRTPGAGGGGWRGGSCTSQR

TRTRRLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:


ACCN: NC_007605

ORF Size: 309 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NC_007605.1](#), [YP_401710](#)

RefSeq ORF: 309 bp

Locus ID: 3783757

UniProt ID: [P03211](#)

MW: 11.1 kDa