

Product datasheet for **VC101201**

UL9 (NC_006273) Virus Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	UL9 (NC_006273) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	UL9
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC101201 represents NCBI reference of YP_081469 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGCAAGAGACTTCAGGTGTTCCCATGGATTACGATCCTTTTCTATACGTCCAAGAGCCAGTATTGGA
ATTATATGACTATTCCTTGACACCTACAGTCGGCTACGGGTCTACAATATCTCACTGCATCCTCTGAA
TAATAGTCTGTTTCAGGACGATGTGTTGAATGGTATATAGACAAACCCATGGTTACTAATAAGCCTTTGT
CTCTACCAGTCCAACGAACGGATCAAGAGCAACCTGGACTCCCCAACATAATGTGGCAGTGTACAGACA
ACCGGACACTGATCCTGATGAATCTGACAACCAGTATTCCAGGAACATTATTTTCAGTCTTTCAAGTA
CCTCGGACAGGGGTGCCAAGCCAATAATCTTTGTTATAACGTCAGTGTCCATTTCAACCCATCAAACC
CACTGTCAACGACCACTTCTCTCTTTACCCTCCAACCTCTGTCCAGACTCTCTCGAAATCAGCCAGA
GTTTTACTAGTACAACTTACACATACC GCCCTCCACTATGCTGCTGGAACGTGGAGGCTCAACATGA
CACTGCAACATCCCATACGATGTGGATTATCCCCTCGTAATAGTGATCACGATTATCGTCTCATCTGT
TTAAGTTCAGAAAGCATGGAATAAGTTTACGCAGTACCGGATTCTGGCATGCTGGCGGCCGC

ACGCGTACGGCGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >VC101201 representing YP_081469
Red=Cloning sites Green=Tags

MSKRLQVFPWITILFYTSKSQYWNMTIPCTPTVGYGSHNISLHPLNNSLFQDDVFEWYIDKPMVTNKLC
 LYQSNERIKSNLDSPNIMWQCTDNRTLILMNLTTTYSRNYFQSFKYLGGQVPKPNNLCYNVSVHFTHQT
 HCHTTTTSSLYPPTSVHDSLEISQSFTSTNFHTAVHYAAGNVEAQHDTATSHMTWIIPLVIVITIIIVLIC
 FKFPQKAWNKFTQYRYSGLAAA

TRRRLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NC_006273

ORF Size: 699 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NC_006273.2, YP_081469
RefSeq ORF: 699 bp
Locus ID: 3077459
MW: 27.0 kDa