

Product datasheet for **VC101211**

UL20 (NC_006273) Virus Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	UL20 (NC_006273) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	UL20
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC101211 represents NCBI reference of YP_081479 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTCGGTATAAGAGCGATGCTTGTAAATGCTCGATTATTATTGGATACAGCTTATCACGAATAACGGCA
CGCGCTCCAATAATACAGACACAATTTTCGTCAGCCTGTTGACCGGCCAAATGGTGTGACACGGACAGC
CATTGGCGGTCTGTACTCTAACTACACAAATTTGACCGGAGCCTTGGATTTACATCTACAACATGAGC
GCAACCAACTCTCCGCTGAAGATAACTGGTCTGTGACAAACCTGACAGAGTCTTGATTAATAAGAGGCG
AAAGCTACGTTACCACTATCTGGTTGCTCGATTGTACAAAAACGATACCTATTGGTACTACGGCAATGC
TTATAACCATACCTGCGAAGGCACGATCAGCGGTATCTGCTTGGCATGTGTAAGTTGTGGAAGAGTTGG
GTGAATAACATCACCTCATACAACACTGTAAGGGTCGAGTCCCTGGGCAACGAAACAAGGTGCATGCTGC
TTCCACGACAGTACACACTGAATGCTACCGTCGAATGGTATAACAAGTCCGAAGGCGATGTCCAGAAGA
GTTTATGGACTACGTATCCTGACTCCCTTGGCAGTACTCACCTGTGGTCTCCAAGAGGCCTATATCCTG
GACAAGGGACGACGCTACATGTACTTGTTCAGCGTGTCTTGCCTGGGGATCACCGGGACTGTGTCTATTA
TTCTGGTTAGCCTCTCCCTGCTGATTCTGATTTGCTACTATAGATGCGGACGATTGTTGATTTGCCCGAG
GGGTTTTGAACCTCTCCCGAGTTCACAGAAGAGGAGGAGGAGAAAGAAAACTTCTCACCCACAACGAT
ATTGAAGTGCAGGTGCCCATACGCACTAGGAGACTGCTGGTGCCTTGGATTAGAGAGAGCAAGATGTGGG
TCTGCCCCACCACTCCCCCGCGGCCCTCACCTCATCGAGTTTCCACCCAGTCCACCCAGCCAGCC
TGGACCTATGCATGGTGGTGTGCATGCCAGCA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >VC101211 representing YP_081479
Red=Cloning sites Green=Tags


MLGIRAMLVMLDYWIQLITNNGTRSNTDTIFVSLLTGPNGVTRTAIGGLYSNYTNLTGAFGFTSTNMS
 ATNSSAEDNWSVTNLTESCINRGESYVTIWLLDCTKNDTYWYYGNAYNHTCEGTISGYLLGMCKLWKS
 VNNITSYNTVRVESLGNETRCMLLPRQYTLNATVEWYNKSEGDVPEEFMDYVILTPLAVLTCGLQEAYIL
 DKGRRYMYLFSVSCVGITGTVSIILVLSLLILICYRCGRLLICPRGFELLPEFTEEEEEKEKLLTHND
 IEVQVPIRTRRLLVPWIRESKMWVLPPLPPRPHLIEFPPSPPPSPGPMHMVVCMPA

TRTRRLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NC_006273

ORF Size: 1014 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NC_006273.2, YP_081479
RefSeq ORF: 1014 bp
Locus ID: 3077434
MW: 38.3 kDa