

Product datasheet for VC101291

UL120 (NC_006273) Virus Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	UL120 (NC_006273) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	UL120
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC101291 represents NCBI reference of YP_081559 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTACCGGGCCGGCGTGACCCTGCTGGTGGTGGCCGTGGTGGAGCTTCGGACGGTGGGACTCCGTGACAG
TCGCCACAACCATTCGGGTCGGCTGGTGGTACGAGCCCCAGGTGAAGATGGCCTACATCTACGAGCACAA
CGACACCAACCTGACAATCTTCTGTAACACCACCGCCTACGACAGCCCTTCTGGCTAGCGGCATGATG
ATCGTGTGCCACCGGACCCAATTCGACCCGGAAGGTGAACACTACGCGAAGACATGGAGAACATCA
AGCAGAAGTACACCCATCAGCTGACCCACATGCTGACCGGCGAGCCCGGCACCTATGTGAACGGCAGCGT
GACTTGCTGGGAAGCAACGGCACCTTCGGCGCCGGAACCTTCATCGTGCGGAGCATGGTGAACAAGACC
GCCGGCAACACCAACACCTTCATCCACTTCGTGGAGGACAGCGAGCTGGTGGAGAACCCCGCCTACTTTC
GGCGGTCCGACCACCGGCCTTCATGATCGTGATCCTGACCCAGGTGGTGTTCGTGGTGTTCATCATCAA
CGCCAGCTTCATCTGGAGCTGGACATTCGGCGGCACAAAAGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:	>VC101291 representing YP_081559 Red=Cloning sites Green=Tags
-------------------	--

MYRAGVTLVAVVVSFGRWDSVTVATTIRVGWVYEPQVKMAYIYEHNDTNLTIFCNTTAYDSPFLASGMM
IVLPHRTQFLTRKVNYSMDMENIKQNYTHQLTHMLTGEPGTYVNGSVTCWGSNGTFGAGTFIVRSMVNK
AGNTNTFIHFVEDSELVENPAYFRSDHRAFMIIVILTQVVFVFIINASF IWSWTFRRHKR


TRTRRLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NC_006273

ORF Size: 603 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NC_006273.2](#), [YP_081559](#)

RefSeq ORF: 603 bp

Locus ID: 3077520

MW: 23.1 kDa