

Product datasheet for **VC101293**

UL122 (NC_006273) Virus Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	UL122 (NC_006273) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	UL122
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>The Viral ORF clone VC101293 represents NCBI reference of YP_081561 with codon optimized for human cell expression

Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGAGAGCAGCGCCAAACGGAAAATGGATCCTGACAATCCCGATGAGGGACCTAGCAGCAAAGTGCCAC
GGCCCGAGACCCAGTGACCAAGGCCACCACATTCCTGCAGACTATGCTGCGGAAGGAGGTCAACAGCCA
ACTGTCCCTGGGCGACCCCTCTCCCCGAGCTCGCCGAAGAGAGCCTGAAAACCTTCGAGCAGGTCAAC
GAGGACTGCAACGAGAACCCGAAAAGGACGTGCTGACCGAGCTCGGCGACATCTCGCACAGGCCGTGA
ACCACGCCGCATCGACAGCAGCAGCACCAGCCACACTGACCACCCACAGCTGCAGCGTGTCCAGCGC
CCCCCTGAACAAGCCACTCCCACCTCCGTGGCCGTACCAATACCCCTCGCCGGCGCCTCCGCCACC
CCCGAGCTGAGCCCCGGAAGAAGCCCCGAAAACCCACACGGCCCTTCAAGGTCATCATCAAGCCCCAG
TGCCACCCGCCCATCATGCTCCCTGATTAGCAGGAGGACATTAAGCCCGAGCCCGACTTCACCAT
CCAGTACCGGAACAAGATCATCGACACTGCCGGCTGCATCGTGATCTCCGACTCCGAGGAGGAGCAGGGC
GAGGAGGTCGAGACCCGGGGAGCCACCGCCAGCAGCCCTCCACCGGCAGCGGCACACCCGGGTGACCA
GCCCCACCCACCTCTGAGCCAGATGAACCACCCCTCCCGACCCCTCGGGCGGCCGATGAGGA
CAGCAGCAGCAGCTCCAGCAGCTGCAGCAGCGCCAGCAGCAGCGAGTCCGAAAGCGAGGAGATGAAG
TGCAGCAGCGGGGGCGCGCAAGCGTGACAAGCTCCACACGGCCGGGGCGGGTTCGGGGCGCCGCCA
GCTCCTCCCTGCTGAGCTGCGGCCATCAGAGCTCCGGCGGCCAGCACTGGCCCCGAAAAAAGAGAG
CAAGCGGATCAGCGAGCTCGACAATGAGAAGGTGCGGAACATCATGAAGGACAAGAATACCCCTTCTGC
ACCCAAAACGTGCAAACACAGCGGGCCGGTGAAGATCGACGAGGTGAGCCGGATGTTCCGGCACACCA
ACCGGTCCTGGAATAACAAGAATCTGCCCTTTCATGATCCCAGCATGCACCAGGTCTCGAGGAGGCCAT
CAAGGTCTGCAAGACCATGCAGGTGAACAACAAGGCATCCAGATCATTTACACCCGGAACCACGAAGTG
AAGAACGAGGTCGACCAGGTGAGTCCCGCTCGGCAGCATGTGCAACCTGGCCCTGAGCACCCCTTCC
TCATGGAGCACACCATGCCCGTACCCACCCCGGACGTGGCACAACGGACCGCCGACGCTGCAACGA
CGGCGTCAAGGCCGTGTGGAACCTCAAGGAGCTCCATACCCACAGCTGTGCCCCGAGCAGCGACTAC
CGGAACATGATCATCCACGCCACACCCGTGGACCTGCTGGCGCCCTGAACCTCTGCCTCCCTCTGA
TGCAGAAGTCCCAAAGCAGGTGATGGTGGGATCTTCAGCACCAACCAAGGGGATTATGCTCCCTAT
TTACGAGACCGCCGCAAAGCCTACGCTGTGGCCAATTCGAAAAGCCACCGAAACCCCGGAAAGAC
CTGGACTCTGTCCCTGGCCATCGAGGCGCCATCCAGGATCTCCGGAATAAGAGCCAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>VC101293 representing YP_081561

Red=Cloning sites Green=Tags

MESSAKRKMDDPNDPDEGPSSKVPREPVTKATTFQLQMLRKEVNSQLSLGDPLFPELAEESLKTFEQVT
EDCNENPEKDVLTGDLAQAVNHAGIDSSSTGPTLTHSCSVSSAPLNKPTPTSVAVTNTPLPGASAT
PELSPRKKPRKTRPFKVIKPPVPPAPIMLPLIKQEDIKPEPDFTIQYRNKIIDTAGCIVISDSEEEQG
EEVETRGATASSPSTGSGTTPRVTSPTPLSQMNHPPPLPDPLGRPDEDSSSSSSSSSSASDSESESEEMK
CSSGGGASVTSSHHRGGFGGAASSLLSCGHQSSGGASTGPRKKKSKRISELDNEKVRNIMKDKNTPFC
TPNVQTRRGRVKIDEVSRMFRHTNRSLEYKNLPMIPSMHQVLEEAIKVCKTMQVNNKGIQIYTRNHEV
KNEVDQVRCRLGSMCNLALSTPFLMEHTMPVTHPPDVAQRTADACNDGVKAVWNLKELHTHQLCPRSSDY
RNMIHAATPVDLLGALNLCPLMQKFPKQVMVRIFFSTNQGGFMLPIYETAAKAYAVGQFEKPTETPPED
LDTLSLAIEAAIQDLRNKSQ

TRTRRLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:


ACCN: NC_006273

ORF Size: 1740 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NC_006273.2, YP_081561](#)

RefSeq ORF: 1740 bp

Locus ID: 3077563

MW: 63.1 kDa