

Product datasheet for **VC101377**

U29 (NC_001664) Virus Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	U29 (NC_001664) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	U29
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC101377 represents NCBI reference of NP_042922 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAACTCCAAGTCTAGTGCTCGCGCCGCATAGTAGATACCGTGGAGGCCGTGAAAAGCGGAAATACA
TCTCTATTGACGAGGGCACACTGAACAATGTCGTGGAGAAGGAGCGAAAGTTCTGAAACAATTTTTGAG
CGGCCGCCAGAACCTGCGAATCGCTGCAAGAGTATTTACCCCTTGTGAGCTGCTGGCACCAGAACTGGAA
AACCTGGGAATGCTGATGTATAGATTTGAGACCGATGTTGACAATCCCAAATCCTGTTTGTGGCCTCT
TCTTCTCTGTAGTAATGCCTTTAACGTGAGTACTTGCCTCGGACAGCCCTCACGGCTATGTATACCAA
CTCCATGGTTCGATAATGTCCTTTCTATGATTAATACCTGTCGCTATCTTGAAGATAAGGTGTCAGTGTT
GGCGTCACAAGCCTCGTGTATGTGGGAGCTCATGTCTTGTCTGTGTATGCAGGGGAATGTGTACG
ATGTTAATAAGGAGAACATCTACGGTCTGACTGTTCTGAAGGAGATAATACTGGAACCCGACTGGGAGCC
AAGGCAGCATTCCACGCAGTATGTGTACGTCGTCATGTATAACAAGGAGGTGCTGGCGAAACTGCAGTAC
GGTATTTACGTGGTCCTTACATCATTCCAGAATGAGGACTTGATCGTCGATATCCTGAGACAGTACTTTG
AAAAGAGCGGTTTCTGTTTTTGAACCTGATCAATTCCAACACCACCCTGAGTACTCGGGTCAGT
TCAGCGCATTGGACGATGTGCCACGAAGACATCAAGAGCGGTTTCTTGCAGTATCGGGGATAACCTTG
TCAGTCATCAAGCTGAAAATATCTTCGTGACCTTTTCAGAGAAGAAAGTTTTTGTCT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >VC101377 representing NP_042922
Red=Cloning sites Green=Tags

MNSKSSARAAIVDTVEAVKKRKYISIDEGTLNNVVEKERKFLKQFLSGRQNLRIARVFTPCCELLAPELE
 NLGMLMYRFETDVPNKILFVGLFFLCSNAFNVSTCVRTALTAMYTNMVDNVLMSINTCRYLEDKVSLE
 GVTSLVSCGSSCLLSCVMQGNVYDVNKENIYGLTVLKEIILEPDWEPRQHSTQYVVVVHVYKEVLAKLQY
 GIYVVLTSFQNEDLIVDILRQYFEKERFLFLNYLINSNTLSYFGSVQRIGRCATEDIKSGFLQYRGITL
 SVIKLENI FVDLSEKKVFV

TRTRRLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NC_001664

ORF Size: 897 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NC_001664.2, NP_042922
RefSeq ORF: 897 bp
Locus ID: 1487906
MW: 34.2 kDa