

## Product datasheet for **VC101439**

### DR6 (NC\_001664) Virus Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	DR6 (NC_001664) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	DR6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC101439 represents NCBI reference of NP_043000 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGAGCACTCGACATACCCAGAGAACAGGGGGACGCATTTCCGTTCCGGCGGAACGAAGCCCGGCTGGCGC  
ACGCCCCAGCCAGAGCCAGGTTTCGATTGGTTGCTGTTGGCACGGGGCCGCCCTCCAAGCTGATGGCTA  
CGCTTTCGCCATCGCGGTGAGCTGATCCACCTTCCCTGGCTCCCTGCTGGTGTCTTGAGCTCCACCAC  
GATCCTTATAGAGATGCTAGATCTGCAACAGTGTGGGGCCACAGATGGGGCTGGCCTGCCACCCAGTGC  
GGCCCCGCTGTGTTTCAGGATTGCGCCCTCGACAGTTGCCTCTACGTGTGTTGTGGGTACGGGAAAAGTT  
GCAGCCTGTCCGTTTCGTGAAGTCATATGTGACGAACAGTCAGCTGGATACGCTTAGAGTGTGCTCGTT  
GGGAAGGACGGTGCCGTCTACGTGCATCACATGCGGGCAGCCGACTTTGCCGCTTGCTAGCAGCACCA  
CAGAGTTCACACGACGGGGCTCCAGCGCGATGCTGTACCTATGAGGAGGATCTTGAACCTCCAGACCA  
GCGGATGTGTGGCACTAATGCACGGCATCTTTGACGTTATCGCGGGCGGCAGACGAACACAACCTG  
CTGACAGTCGGGGCCTCTGCCAGACACATGCAGGAGTAAGCTGCAATCTGCTCGAAACCGTCGGGGATC  
CATGGACAGCCGTGCCTGCCGAAGGATGACTTTGACCGTGCCTCAGGTGCAGTACCGATTGTGGCCGA  
GGCCCCAAGAGATCTGCGGAGGCATCTTATGCAGGGCACCCACTCGGTCTTGGTTGGTCTGCGGCGTG  
CTGTCCCAGGAGAGAACTCAAAAACCGAGCCACCGATCAGGACTACTGTCCGGAACGTGCCCACTC  
CCGGTCTCGGGAAGTTGAAATCGCTTGGGTGGTGTGACTTTGGCCGGGCCACTTCTGGCATTGTTGGCC  
TGATACAGGTAAGATATGCAGGCTCGCAAATTCATTGACACCCTGTGGAAAATGGGACCTCGGGCAATG  
AGAGGACATTGGACTTATAGTGCACCTGGCCGCCACCTGCCGGGTGATGCATGGCCCTTGTGTGAGCAGC  
TGAGGCCTCAGTGAAGGCATCTGAAGAAGAGTCTGTGCCTGGACTGGACGCTACGGAAGTGAAG

**ACGCGT**ACGCGGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


**Protein Sequence:** >VC101439 representing NP\_043000  
Red=Cloning sites Green=Tags

MSTRHTQRTGGRISVRRNEARLAHARARARFDWLLLAGRPSKLYGYASRHRGELIHLPPWPCWCLELHH  
 DPYRDARSATVWHRWGWPATHVRPRCVDCLDSCLYVCCGYGEKLQPVGFVKSYYVTNSQLDTRLRVLLV  
 GKDGAVYVHHMRAARLCRLASSTTEFTRRGLQRDAVTYEEDLELPDQRMCGTNARHLFDVIAAAADEHNL  
 LTVGGLCQTHAGVSCNLLLETVGDPWTAVPAARMTLTVPQVQYRLWPEARDDLRRHLYAGHPLGPWLVCV  
 LSRERETQKPSPIRTTVGNVPTPGPREVEIAWVVLTLAGPLLAFWPDTGKICRLANSFSTLWKMGP  
 RGHWTYSAPGRHLPDGAWPLCEHVRPQWKASEEESVPGDATEVK

TRTRRLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NC\_001664

**ORF Size:** 1185 bp

**OTI Disclaimer:** The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NC\\_001664.2](#), [NP\\_043000](#)

**RefSeq ORF:** 1185 bp

**Locus ID:** 1487974

**MW:** 44.6 kDa