

Product datasheet for **VC101615**

U75 (NC_001716) Virus Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	U75 (NC_001716) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	U75
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC101615 represents NCBI reference of YP_073815 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGACTGCCCCGAAGCCTGAAGGGATTTGAGACTCTGGATATTTCAACATCAGTTCTGTTAGAAACA
TTTGTGAGCTGAATAAGCATGCCGAACACGTCTCAATTAATAATTTGCCCTTGATCGACATCAGTGTGG
CAACGACGACGTGTGGTTCCACCTTGAAGATGGAACATCGTAAATGGAAAGACTATAAAAGTATCTGT
GAGAAAACGCTGGTTTTCTGGCTTCATTGGGATCATCTCCTGGATAGCGAAGACACTTTGAGGAAA
TTCGGCTGTCAAAGACCCAGTGTAAACGCGCATTATTTATTTGATTCTCAAGGAAGATACAGAGTTCCT
GCTTTGCGGAATCGTCTACGCTTTGGAGAATCTGCCCATCAAGGGGCAGACTCTCCTGAAGCTTCGGGAC
ATTATCAAGAAGATTAGTGTGACTTCTGTTAGCAAGCTTCTCGCTTGACATGTCAGAAGCTCATT
CCATCCTGCGATATCTTTTACGACGACAAACAGCAGGACGTCCTGGATAAGGTCCCACCAATCATTCA
GTTGTAACGAGTCTAAAAGTCTAATATCCATATGCTGAACCTGTTCTTCAAATCTCAGGACAACGAC
GACTTGTACTCTGTCCCTGAATACCAGACGGCTCCAGGACGATTCTAAGTATCTGATCGATTTTCTTA
AGAGTGCCATTTGCGACGCTTTCAGTAAAGAGTATAAAATGACGGAATCGAAAAAATCTCCCTCCAC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >VC101615 representing YP_073815
 Red=Cloning sites Green=Tags

MTLPASLKGFE~~LD~~IFNISSVRNICELNKHAEHVS~~IKNLPLIDISVGNDDVWFHLEDGTIVNGKSYKSIC~~
 E~~KT~~L~~G~~F~~L~~G~~F~~I~~G~~I~~I~~L~~L~~D~~S~~E~~D~~T~~L~~E~~E~~I~~R~~L~~S~~K~~T~~Q~~C~~K~~R~~R~~I~~I~~Y~~L~~I~~L~~K~~E~~D~~E~~F~~L~~L~~C~~G~~I~~V~~Y~~A~~E~~N~~L~~P~~I~~K~~G~~O~~T~~L~~L~~K~~L~~R~~D
 I~~I~~K~~K~~I~~S~~V~~T~~L~~P~~V~~S~~K~~L~~L~~A~~C~~T~~C~~Q~~K~~L~~I~~S~~I~~L~~R~~Y~~I~~F~~Y~~D~~D~~K~~Q~~Q~~D~~V~~L~~D~~K~~V~~P~~P~~I~~I~~Q~~L~~Y~~E~~S~~K~~T~~A~~N~~I~~H~~M~~L~~N~~L~~F~~F~~K~~S~~H~~D~~N~~D
 D~~T~~C~~T~~L~~S~~L~~N~~T~~R~~R~~L~~Q~~D~~S~~K~~Y~~L~~I~~D~~F~~L~~K~~S~~A~~I~~C~~D~~A~~F~~S~~K~~E~~Y~~K~~M~~T~~E~~I~~E~~K~~T~~S~~L~~H

TRTRRLE~~Q~~K~~L~~I~~S~~E~~E~~D~~L~~A~~A~~N~~D~~I~~L~~D~~Y~~K~~D~~D~~D~~D~~D~~K~~V~~

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NC_001716

ORF Size: 768 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NC_001716.2, YP_073815](#)
RefSeq ORF: 768 bp
Locus ID: 3289533
MW: 29.4 kDa