

Product datasheet for VC101625

U91 (NC_001716) Virus Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: U91 (NC_001716) Virus Tagged ORF Clone
Tag: Myc-DDK
Symbol: U91
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >The Viral ORF clone VC101625 represents NCBI reference of YP_073825 with codon optimized for human cell expression
 Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTACACCCCTGGAGTATGAGAAGAGGGTATCCCGCCCTAAATTGACATATTGGATCATCCTGGCCATCT
 TGTTTCGTCTTTCTATCATCACAGGCTCCGTTCTGATCGTGATAGAGACCCTTTCTATTCAGAGGACAAC
 CCTGAATGCACAAAACGACAAGACATCTACAGTGGTGCCGGAACCTACCAGCAACAGTCCCGATCAAACC
 ACCGTGACAAACTTTTCAGCGAGCAGCAAGCCAACCTTTGTCTAGCAAGCAGCCCGTTGGATTAGGCC
 TTACCACAGCTTTCGGCATCCTCACTCTGTTTAGTGTTATGATGATCATTATTACTTGAATTTTTGGCT
 GACAGAAAAGAACGACAAGACCGCCAACCTACCGAATATTATTCTGAAGACATCCTCGACTACACAAAC
 CCCAGTTTTACAGAAATCGACGAAGATTCTCTAAGGTG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >VC101625 representing YP_073825
 Red=Cloning sites Green=Tags
 MYTLEYEKRVSRPKLTYWIIILAILFVFLIITGSVLIVIETLSIQRTTLNAQNDKTSTVVPELTSNSPDQT
 TVTNFSASSKPTLSSKQPGWIQALTTAFGILTLFSVMMIITCNFWLTEKNDKTANPTEYYSEDILDYTN
 PSFTEIDEDSSKV

TRTRRLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:


ACCN: NC_001716

ORF Size: 459 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NC_001716.2](#), [YP_073825](#)

RefSeq ORF: 459 bp

Locus ID: 3289543

MW: 17.3 kDa