

Product datasheet for **VC101630**

K1 (NC_009333) Virus Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	K1 (NC_009333) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	K1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC101630 represents NCBI reference of YP_001129350 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGTTTTGTATGTGGTTTGCAGCCTCGCGTGTGTTTTAGAGGACTGCTTCTCTGTCAGTCTA
GCCCCAACCTCTGCCAGGCGTGATCAGTACCCCATACACGCTCACATGTCCCTCAACACAAGCCTGCC
CACTTCTTGGTACTGTAATGACACTCGTTGTTGCGAGTGACACAGGGCACATTGACAGTGGACACATTG
ATTTGCAATTTCTCATGCGTAGGACAGAGCGGGCACCCTACAGTCTTTGGATTACATGGTATGCCAAC
CGTTCTTCAGACCTTTTGCGGCAGCCAGTAATACGGTGACATGCGGCCAGCACGTCACCTTTATTG
CTCTACTTCCGGTAACAATGTGACAGTGTGGCATCTGCCAAACGGCCAGAACGAAACCGTCTCACAGACT
AAATATTACAACCTCACTCTGATGAACCAACCGAAGGGTCTATGCATGCAGTAACGGGCTTTCTCAA
GGCTGTCTAACCGGCTGTGCTTTTCTGCTAGATGCGCCAATCACTCCAGAGACACATACAGTCTCTGT
GAGCAGCACCACGGGATTTGGACTTTTGGCACTGCCCCACACTGTTCTGTTATGAAAGAGGTCAAATCA
ACCTATCTGTACATTCAAGAGCATCTTCTCGTGTATGACACTCGTCCGACTGATTGGCACCATGTGCG
GGATCCTTGAACCATTATCTTTGCCACTGTCAAAGCAGAGGGACAGTAACAAGACTGTTCCCAACA
GCTGCAGGATTACTATAGCCTGCATGATCTTTGTACGGAAGATTACACTCAGCCAGTCTGACTGGTAC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence: >VC101630 representing YP_001129350
 Red=Cloning sites Green=Tags

MFLYVVVCSLAVCFRGLLSLSLQSSPNLCPGVISTPYTLTCPSTSLPTSWYCNDTRLLRVTQGLTVDTL
 ICNFCVCGQSGHRYSLWITWYAQPVLQTFCGQPSNTVTCGQHVTLYCSTSGNNVTVWHLPNGQNETVSQT
 KYYNFTLMNQTEGACYACSNGLSSRLSNRLCFSARCANITPETHTVSVSSTTGFRFTATAPTLFVMKEVKV
 TYLYIQEHLVFMFLVALIGTMCGLGTIIFAHQCQQRDSNKTPQQLQDYYSLHDLCTEDYTPVDWY

TRTRRLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NC_009333

ORF Size: 837 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NC_009333.1](#), [YP_001129350](#)
RefSeq ORF: 837 bp
Locus ID: 4961511
UniProt ID: [Q2HR82](#)
MW: 31.2 kDa