

Product datasheet for **VC101713**

ORF74 (NC_009333) Virus Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ORF74 (NC_009333) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ORF74
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC101713 represents NCBI reference of YP_001129433 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTGCCGAAGATTTTCTGACAATATCCTGGACGACGATGAAAGCTGGAACGAGACCCTCAATATGA
GCGGCTATGACTATCCGGGAATTTTCCCTCGAAGTGAGCGTGTGCGAGATGACCACTGTGGTGCCTTA
CACTTGGAAACGTTGGGATCCTGAGCCTCATCTTCTGATCAATGTGTTGGGCAACGGTTTGGTGACCTAT
ATCTTTTGTAAAGCACAGTCCCGCGCAGGCGGATTGACATTCTGCTCCTCGGCATCTGTTGAATTCTC
TGTGCTTGTCCATTTCACTGCTGGCCGAGGCTTGATGTTTCTTTCCCTAATATCATCAGCACTGGCCT
GTGACGGTTGGAGATCTTCTTTTATTACTTGATGTTTATTTGGATATCTTCAGCGTCGTCTGCGTGTCT
CTGGTCAGGTACCTGTTGGTTGCTTATAGTACACGAAGCTGGCCTAAAAAGCAAAGTCTTGGATGGGTGT
TGACCAAGTGTGCTTTGTTGATCGCTCTGGTCTGTGTCAGGAGACGCTGTAGACACCGAAGCAGAGTTGT
TGACCCGGTGAGTAAGCAGGCCATGTGTTATGAGAATGCAGGAAACATGACCGCTGACTGGAGGTTGCAT
GTCCGGACAGTGTGAGTACGGCGGGTTTTCTCCTGCCTCTCGCTCTTCTGATCCTTTTTATGCACTGA
CCTGGTGCCTTGTTCGAAGAACTAACTTCAGGCCAGACGGAAGTTAGGGGAGTGATAGTGGCTGTGGT
GCTGTTGTTTTTGTGTTCTGCTCCCGTACCAGTGCTGAACCTGCTTGATACATTGCTTCGCCGCGG
TGGATCCGGGATAGTTGCTACACCAGAGGACTGATTAACGTGGGACTGGCGGTAACCTCACTTTTGCAGG
CACTCTACTCCGCTGTCGTGCCTCTGATTTACAGTTGTTTGGGCTCTCTCTCCGGCAGAGGATGTACGG
TCTGTTTCAAAGCCTGCGCCAGTCTTTCATGAGCGGAGCAACCACC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >VC101713 representing YP_001129433
Red=Cloning sites Green=Tags

MAAEDFLTIFLDDDESUNETLNMSGYDYSGNFSLEVSVCEMTTVPYPTWNVGILSLIFLINVLGNGLVTY
 IFCKHRSRAGAIIDILLGICLNSLCLSIALLAEVLMFLFPNIISTGLCRLEIFFYYLYVYLDIFSVCVS
 LVRVLLVAYSTRSWPKQSLGWLTSAALLIALVLSGDACRHRSRVDPVSKQAMCYENAGNMTADWRLH
 VRTVSVTAGFLLPLALLILFYALTWCVRRRTKLQARRKVRGVIVAVLLFFVFVCFYPYHVLNLLDTLRRR
 WIRDSCYTRGLINVGLAVTSLQLALYSVVPLIYSCLGSLFRQMYGLFQSLRQSFMSGATT


TRTRRLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/ja2476_b04.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NC_009333

ORF Size: 1026 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NC_009333.1</u> , <u>YP_001129433</u>
RefSeq ORF:	1026 bp
Locus ID:	4961460
UniProt ID:	<u>Q2HR82</u>
MW:	38.7 kDa