

## Product datasheet for **VC101717**

### rev (NC\_001802) Virus Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** rev (NC\_001802) Virus Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** rev  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >The Viral ORF clone VC101717 represents NCBI reference of NP\_057854 with codon optimized for human cell expression  
**Red**=Cloning site **Blue**=ORF **Green**=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCAGGGAGAAGTGGCGACAGTGATGAAGAACTGATCCGGACCGTCCGGCTGATTAAGCTCCTGTATC  
AGAGCAACCCCCACAAATCCAGAGGGGACCAGACAGGCAAGGAGGAATAGGAGACGCAGATGGCGGGA  
GCGACAGCGCCAGATCCATAGTATCAGTGAGAGAATCCTCGGCACTTATCTCGGGAGATCAGCAGAGCCC  
GTCCCTCTCCAGCTCCCACCTCTTGAACGATTGACCCTGGACTGTAACGAGGATTGTGGTACGAGCGGAA  
CCCAGGGCGTCGGCAGTCCACAGATTCTCGTGGAAAGCCCGACTGTGCTCGAATCAGGGACAAAGGAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >VC101717 representing NP\_057854  
**Red**=Cloning sites **Green**=Tags  
MAGRSGDSDEELIRTVRLIKLLYQSNPPNPEGTRQARRNRRRRWRERQRQIHSISERILGTYLGRSAEP  
VPLQLPPLERLTLDNEDCGTSGTQGVGSPQILVESPTVLESATKE  
**TR**TRRLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mk8100\\_a02.zip](https://cdn.origene.com/chromatograms/mk8100_a02.zip)

**Restriction Sites:** Sgfl-MluI



[View online »](#)

Cloning Scheme:



ACCN: NC\_001802

ORF Size: 348 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NC\\_001802.1](#), [NP\\_057854](#)

RefSeq ORF: 348 bp

Locus ID: 155908

UniProt ID: [P04585](#)

MW: 13.1 kDa