

## Product datasheet for **VC101755**

### M (NC\_004148) Virus Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	M (NC_004148) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	M Protein
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC101755 represents NCBI reference of YP_012607 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGC**C

ATGGAGTCATACCTCGTGGACACGTACCAGGGGATACCTTATACCGCCCGCTGCAAGTAGACCTCGTGG  
AAAAAGACCTCCTGCCCGCCTCCCTGACTATCTGGTCCCCTCTTTAGGCAAATACACCTCCTGCAGT  
GCTGCTGGATCAACTGAAGACACTCACTATAACCACCCCTCTATGCCGCCAGCCAAAGTGGACCTATTCTG  
AAAGTAAACGCATCAGCTCAAGGGCCGCCATGTCTGTGCTTCCTAAGAAATTTGAGGTCAACGCGACCG  
TCGCTCTGGATGAATACTCCAAGCTGGAGTTCGATAAGCTGACGGTCTGCGAAGTGAAGACTGTGTACCT  
TACCACTATGAAGCCATATGGCATGGTGGAGCAAATTCGTCTCAAGCGCTAAACCCGTGGCAAAAAAACA  
CATGATCTGATCGCCCTGTGCGATTTTCATGGACTTGAAAAGAATACACCTGTGACCATTCCGGCCTTCA  
TTAAGTCCGTGAGCATTAAAGGAATCTGAGTCCGCGACAGTGAAGCAGCAATCAGTTCAGGAGCCGACCA  
GGCCCTTACCCAAGCGAAAATTGCCCTTATGCAGGACTTATCATGATCATGACCATGAATAATCCCAA  
GGCATTTCAAAAAACTGGGCGCTGGAACACAGGTGATTGTAGAGCTCGGCGCTTACGTCCAGGCCGAGT  
CTATAAGCAAGATCTGTAAAACCTGGTACACACAGGCCACACGGTACGTTCTCAAAGCAGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


**Protein Sequence:** >VC101755 representing YP\_012607  
 Red=Cloning sites Green=Tags

MESYLVDTYQGIPYTAAVQVDLVEKDLLPASLTIWFPLFQANTPPAVLLDQLKTLTITTLAASQSGPIL  
 KVNASAQGAAMSVLPKKFEVNATVALDEYSKLEFDKLTVCVKTVYLTTMKPYGMVSKFVSSAKPVGKKT  
 HDLIALCDFMDLEKNTPTIPAFIKSVSIKESSEATVEAAISSEADQALTQAKIAPYAGLIMIMTMNNPK  
 GIFKKLGAGTQVIVELGAYVQAESISKICKTWSHQTRYVLKSR

TRTRRLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NC\_004148

**ORF Size:** 762 bp

**OTI Disclaimer:** The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NC\\_004148.2](#), [YP\\_012607](#)  
RefSeq ORF: 762 bp  
Locus ID: 2799938  
MW: 27.6 kDa