

## Product datasheet for VC101925

### hypothetical protein (NC\_001457) Virus Tagged ORF Clone

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Expression Plasmids  |
| Product Name:             | hypothetical protein (NC_001457) Virus Tagged ORF Clone  |
| Tag:                      | Myc-DDK  |
| Symbol:                   | hypothetical protein   |
| Mammalian Cell Selection: | Neomycin   |
| Vector:                   | pCMV6-Entry (PS100001)   |
| E. coli Selection:        | Kanamycin (25 ug/mL)   |
| ORF Nucleotide Sequence:  | >The Viral ORF clone VC101925 represents NCBI reference of NP_040889 with codon optimized for human cell expression<br>Red=Cloning site Blue=ORF Green=Tags(s) |

GACGTTGATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCCGACGGACGCCGACTCTGGACGACTTTTGCAGGCGCTTCGACATTAGCTTTTTCGATTTGC  
GCCTTACTTGTATCTTCTGCAGTCACACCGTTGACCTGGCGGATTTGGCTCTGTTCTACTGAAGAAGCT  
TTCCCTGGTTTTCCGGGCAACTGCTACTATGCCTGCTGCTCCGAATGCCTCAGACTTCCGCTCTCTTT  
GAACAAGAAAATTACTTCCAGTGTCCATTAAGGCCGTCACCTCGAGGAGATTGCACAGAAAAAATTA  
AAGAGATTTGCATAAGATGTATCTGTTGCCTGCGCTGCTTGATATTGTGGAGAAGTTGGATCTCTGT  
CTCAGACGAAACCTGCTATCTGATCCGGGGCTCTGGCGGTTATTGCCGAACTGTATACGCAAGCAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

|                   |  |
|-------------------|--|
| Protein Sequence: | >VC101925 representing NP_040889<br>Red=Cloning sites Green=Tags |
|-------------------|--|

MADGRPATLDDFCRRFDISFFDLRLTCIFCSHTVDLADLALFYLLKLSLVFRGNCYYACCSECLRLSALF  
EQENYFQCSIKAVHLEEIAQKKIKEICIRCICLLRLLDIVKLDLLYSDETCYLIRGLWRGYCRNCIRKQ

**TR**TRRLEQKLISEEDLAANDILDYKDDDDKV

|                    |           |
|--------------------|-----------|
| Restriction Sites: | Sgfl-MluI |
|--------------------|-----------|



[View online »](#)

## Cloning Scheme:



ACCN: NC\_001457

ORF Size: 420 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NC\\_001457.1](#), [NP\\_040889](#)

RefSeq ORF: 420 bp

Locus ID: 1489450

MW: 16.5 kDa