

## Product datasheet for **VC101932**

### ORF E6 (NC\_001354) Virus Tagged ORF Clone

#### Product data:

|                           |                                                                                                                                                                              |
|---------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:             | Expression Plasmids                                                                                                                                                          |
| Product Name:             | ORF E6 (NC_001354) Virus Tagged ORF Clone                                                                                                                                    |
| Tag:                      | Myc-DDK                                                                                                                                                                      |
| Symbol:                   | ORF E6                                                                                                                                                                       |
| Mammalian Cell Selection: | Neomycin                                                                                                                                                                     |
| Vector:                   | pCMV6-Entry (PS100001)                                                                                                                                                       |
| E. coli Selection:        | Kanamycin (25 ug/mL)                                                                                                                                                         |
| ORF Nucleotide Sequence:  | <p>&gt;The Viral ORF clone VC101932 represents NCBI reference of NP_040285 with codon optimized for human cell expression</p> <p>Red=Cloning site Blue=ORF Green=Tags(s)</p> |

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGCCTCAACGTCAGGGGTCGGCTCCGTGGGACCCGCTAGCTGCTGTGAAACCCAAAAACCCACACGA  
 TTCGAGAACTGTGCCTCGCACAGCAGATTACTTATCCATGCATTCAATTGTGCTGCTACTACTGTATAA  
 AATTCTCTCCGTACTGGACATTTATGCCTTTGACCAGAGCTGTCTGTATCTGTATGGGGGGAGGGCGGA  
 CCCACAGGGATTGTTCACAGTGCACCCGCTGCTGGCAAGACTCGAATTTACGGCTAGGCATGAAGTGT  
 CCTGTGCTGCCTCTCGGCTGCCTCACTTTATCGGGCAGTCACTGTCAGATCTCGAAGTTAGGTGTGTCAG  
 ATGCCTTGCTTTGCTGCAGAGTGTGAGAAAGATTATATACTCCGCGAAGATCTGTCCTGCACCGGATC  
 GGGGGGATCTGGCGGGGAACCTGCGTGCCTGCATGGTAGGGTTGTAT

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA


|                   |                                                                                |
|-------------------|--------------------------------------------------------------------------------|
| Protein Sequence: | <p>&gt;VC101932 representing NP_040285</p> <p>Red=Cloning sites Green=Tags</p> |
|-------------------|--------------------------------------------------------------------------------|

MASTSGVGSVGPASCCETQKPHTIRELCLAQQITYPCIQLCCHYCYKILSVLDIYAFDQSCLYLSWEGGG  
 PTGICSQCTRVRLARLEFTARHEVSCASRLPHFIGQSLSDLEVRVRCALLQSVEKDYILREDLSVHRI  
 GGIWRGTCVRCMVGLY

**TR**TRRLE**QKL**ISEEDLAANDILDYKDDDDKV

|                    |           |
|--------------------|-----------|
| Restriction Sites: | SgfI-MluI |
|--------------------|-----------|


[View online »](#)

**Cloning Scheme:**


**ACCN:** NC\_001354

**ORF Size:** 468 bp

**OTI Disclaimer:** The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NC\\_001354.1](#), [NP\\_040285](#)

**RefSeq ORF:** 468 bp

**Locus ID:** 1489281

**MW:** 17.3 kDa