

Product datasheet for **VC101968**

E4 (NC_001691) Virus Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: E4 (NC_001691) Virus Tagged ORF Clone
Tag: Myc-DDK
Symbol: E4
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >The Viral ORF clone VC101968 represents NCBI reference of NP_043427 with codon optimized for human cell expression
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

CTGGTGATACAATATACGCTTAACAGCTTTAACCTCATGCCTTTGTACATGGCCAATCTGGACAATGGGC
TGCTGTACCTGAAGACAAAACCTGTTTATCCCCTGAGCTTGGCTCCCCAGGGGTCCGGCTTGGGCTGCT
CACCAACTCTCCCGACTGCCTCACAGGACCACGAGGGCGAACCGAAGGATCTCGAAGCCGTTAACCAG
AAGCCATATCGGACCCCTAACCATCCTCCTCGACACCAGCAGTACGACTTCGACGAGGACGACGAGAAAG
AGAATACCATTCCAACAGATACAGAGTCTCACAATCAAAATTGGGAACCAACCCTGCGACAATTGCTTAA
GAAATGGGAAGAGGACCTGACCCAATTGCAACGAATGGTGTGTCAGGATTTGGATGATTATAAGCGGAAA
CTCGGGATTCACCACTCC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >VC101968 representing NP_043427
Red=Cloning sites **Green**=Tags

LVIQYTLNSFNLMPLYMANLDNGLLYLTKLFIPLSLAPQGVRLGLLTNSPGLPHRTTRANRKDLEAVNQ
KPYRTPNHPPRHQQYDFDEDEKENTIPTDTESHQNWEPTLRQLLKKWEEDLTQLQRMVCQDLDDYKRK
LGIHQ

TRTRRLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



Cloning Scheme:


ACCN: NC_001691

ORF Size: 438 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NC_001691.1](#), [NP_043427](#)

RefSeq ORF: 438 bp

Locus ID: 1403633

MW: 17.2 kDa