

Product datasheet for **VC102011**

E4 protein (NC_010329) Virus Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: E4 protein (NC_010329) Virus Tagged ORF Clone
Tag: Myc-DDK
Symbol: E4 protein
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >The Viral ORF clone VC102011 represents NCBI reference of YP_001672012 with codon optimized for human cell expression
 Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGATAGAAAAACAGCCTTTCCAACCCATAATGGCAACATATTCATTATCCGCACCCGGACCTCTAGTG
 GCATCAAATGTCAGGTAGTTCTGATCACCATGGGGTGTACCATGATGAGGAAGATCCCTATGAAAGGGTA
 TATATCTTTCTATCTGTCACTCTCTTTGAAGAACATGGCAGTTGTGGACAATGGTCTGTGCATCATCAAG
 ATTCTGCTGTTTTGCCCTCCACCAGCCCTGCTCGGGGTCCACCCCTACCCCAAAGCCAATAT
 CCACACCCAGCCTCCTGCAACCCGCGCCCACTCCCTCGGAGAGAGGTCCAGGGGAGCGATTATAAACG
 ACGGTCCGCCAGCCACCTCCAAGAAGGCCACCATATTCGATCTGCAGGAAGAGGATGAAGAGAAAGAG
 AACGAACCGCCGAGCAGCGGAAGAAGAGCAGCCTCCCACTGAAACACTTTCGAGAACCTGTCCAGGT
 TGCTGAGAAGATGGGATCAAGACATCGAGCATCTGAAAGATCTGGTGAGCAGAGACTTGGACGACTTCAA
 GCGCACACTCGGAATACCTCAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >VC102011 representing YP_001672012
 Red=Cloning sites Green=Tags

MIEKTAFPTHNGNIFIIRTRTSSGIKQVVLITMGCTMMRKIPMKGYISFYLSLSLKNMAVVDNGLCIK
 ILLFLPLPALLGGPHSPKPISTPQPPATPRPLPRREVQGSYKRRSAQPPRRPTIFDLQEEDKEE
 NEPPEQREEEQPPTETLSENLRLRRWDQDIEHLKDLVSRDLDDFKRTLGIPIH

TRTRRLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NC_010329

ORF Size: 582 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NC_010329.1](#), [YP_001672012](#)

RefSeq ORF: 582 bp

Locus ID: 5896662

MW: 22.2 kDa