

## Product datasheet for **VC102038**

### **E4 (NC\_005134) Virus Tagged ORF Clone**

#### Product data:

Product Type:	Expression Plasmids
Product Name:	E4 (NC_005134) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	E4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC102038 represents NCBI reference of NP_932323 with codon optimized for human cell expression <b>Red</b> =Cloning site <b>Blue</b> =ORF <b>Green</b> =Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGC**C

ATGCTGGCTGACCTCGAATTGGTTGAAAACGGCAAGTGTGTCTTGACAGGCAAAGTCTCTGCTGCTGC  
TCCCTGCACCACGACGACCTCATACTAAACGGTCCGAAAATGGACTCAGCCCTCCAAGTCTCCGAGGTC  
CGAGCCTCCACTTACTCCGGCCCCGATCGCTTGGGGATCTGGAGCCAAAACCCCAAAGGGGAGGGG  
ACGGACGAGAAAAATTTGGTCTGCAACTGCCTCCATCAGGATCTGGGAAGAAATATCACGATTTGCTCC  
ACGGTAAGCCAAAGAACGATGGACCTGGCGTGCACAGGCAGGGGAGTCCAAGCCCGTCCATGGCGATGG  
GGATAACGGGGAACAGGTCCGGTCCCAGAACCTGAGGCGGGTCTGGCGCAGGAGCAGGTGCCGACGT  
GAACCTGAGGCCGAGGCCAGGGAAACACACAGAACCAGGGGAGGGGAGGTTGAAGGCAGGACCCCGC  
AGGACCCTGAGGTAGACCTCAAAGCGCCGATCTGCTGACCAGCGTAGCATTTTCAGCTCCATAAGTGGGA  
GGGACAGTTTAACAGCTGGTAGAGACTATACTCCGCGACCTGAGGGATTACTGGATGAAGTTGTCAACC  
CCCCAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >VC102038 representing NP\_932323  
 Red=Cloning sites Green=Tags

MLADLELVENGKCVLTGKLSLLLLPAPRRPHTKRSENGLSPTSPRSEPPLTPAPYRLGDLEPKPKGEG  
 TDEKNLVLQLPPSGSGKKYHDL LHGKPKNDGPGVHRQGQSKPVHGDGDNGEPGPVPEPEAGAGAGAGAAA  
 EPEAGGQNTQNRGEVEGRTPQDPEVDLQSADLLTSAFQLHKWEGQFNQLVETILRDLRDYWMKLS  
 T PQ

TRTRRLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NC\_005134

**ORF Size:** 636 bp

**OTI Disclaimer:** The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NC\\_005134.2](#), [NP\\_932323](#)  
RefSeq ORF: 636 bp  
Locus ID: 2943335  
MW: 22.6 kDa