

Product datasheet for **VC102046**

Y1 (NC_003461) Virus Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Y1 (NC_003461) Virus Tagged ORF Clone
Tag: Myc-DDK
Symbol: Y1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >The Viral ORF clone VC102046 represents NCBI reference of NP_604437 with codon optimized for human cell expression
 Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGAGCTCCGACAGCCTTACCTCATCTTATCCTACTAGCCCTCAGAAATTGGAGAAGACCGAGGCCGGTT
 CTATGGTCAGTCCACTACTCAGAAGAAGACCTCCCATCATGCTAAGCCAACGATTACCTAAGACCGA
 ACAAGCCAGAGGAGGCCCAAAATCATCGACCAAGTCCGGGGGTAGAATCCCTCGGAGAGCAGGTCAGC
 CAGAAGCAGAGGCATATGCTTGAGTCTCTGATTAACAAGGTCTACACCGGCCACTGGGAGAGAACTCG
 TGCAGACACTTTTCTGAGAATTTGGGCTATGAAGGAAACCCCGAGTCAACAAAATTCTTCAGATGCG
 GGAAGATATAAGAGACCAGTACTTGCGCATGAAAACCGAGAGATGGCTGCGGACCCTTATCCGGGCAAA
 AAGACTAAATTGCGGATTTCCAAAACGATATGAGGAGGTTACCCCTACCTGATGATGGAGAGAGTGG
 AGCAGATCATCATGGAAGAGGCCTGAAACTCGCAGCCACATTGTGCAGGAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >VC102046 representing NP_604437
 Red=Cloning sites Green=Tags

MSSDSL TSSYPTSPQKLEKTEAGSMVSSTTQKKTSHHAKPTITTKTEQSQRPKIIDQVRGVESLGEQVS
 QKQRHMLES LNKVYTGPLGEELVQTLYLRIWAMKETPESTKILQMREDIRDQYLRMKTERWLR LIRGK
 KTKLRDFQKRYEEVHPYLMMERVEQIIMEEAWKLAHIVQE

TRTRRLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



Cloning Scheme:


ACCN: NC_003461

ORF Size: 543 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NC_003461.1](#), [NP_604437](#)

RefSeq ORF: 543 bp

Locus ID: 935262

MW: 21.3 kDa