

Product datasheet for **VC102076**

proteinase 3C (picornain 3C) (NC_001897) Virus Tagged ORF Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | proteinase 3C (picornain 3C) (NC_001897) Virus Tagged ORF Clone |
| Tag: | Myc-DDK |
| Symbol: | proteinase 3C (picornain 3C) |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| ORF Nucleotide Sequence: | >The Viral ORF clone VC102076 represents NCBI reference of NP_740736 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s) |

GACGTTGATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTCCATATGATGGACAGCTCGAGCATATTATCTCACAGATGGCCTATATTACAGGATCCACAACCG
GTCACATCACTCATTGCGCCGGTTACCAGCAGCAGCAGATTATTCTCCACGGGCACTCTAAAAATATCT
TGAACAGGAGGAAGAGCTGACATTGCACTATAAAAACAAAGTATTCCAATAGAGCAACCTTCAGTCACT
CAAGTGACCTTGGGGGGAAACCTATGGACCTCGCCATCGTCAAGTCAAATTCGCTTTTCGGTTCAAAA
AGAATTCTAAATATTACACCAACAAAATTGGGACAGAGAGATGCTGATCTGGATGACCGAACAGGGCAT
AATTACAAAAGAGGTGCAGAGGGTCCACCATAGCGGTGGGATCAAACCGGGGAAGGGACCGAGAGCACC
AAGACCATCTCTACACTGTCAAGTCTGCAAGGGCATGTGTGGAGGTCTGCTCATTAGCAAAGTGAAG
GAAATTTCAAGATTCTGGGTATGCATATAGCCGGAATGGTGAGATGGGTGATAGCTATCCCATCAACT
CCTGAAGAATGATATGAGTGATCAG

ACGGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

| | |
|-------------------|--|
| Protein Sequence: | >VC102076 representing NP_740736 Red=Cloning sites Green=Tags |
|-------------------|--|

MAPYDQGLEHIISQMAYITGSTTGHITHCAGYQHDEIILHGHSIKYLEQEEELTLHYKNKVFPIEQPSVT
QVTLGGKPMDLAIVKCKLPFRFKKNSKYTNKIGTESMLIWMTEQGIITKEVQRVHHSGGIKTREGTEST
KTI SYTVK SCKGMCGLLISKVEGNFKILGMHIAGNGEMGVAIPFNFLKNDMSDQ

TRRRLEQKLISEEDLAANDILDYKDDDDKV




[View online »](#)

Chromatograms: https://cdn.origene.com/chromatograms/ja2024_e01.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NC_001897

ORF Size: 585 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NC_001897.1](#), [NP_740736](#)

RefSeq ORF: 585 bp
MW: 21.8 kDa