

Product datasheet for **VC102080**

NS (NC_000883) Virus Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NS (NC_000883) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	NS
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>The Viral ORF clone VC102080 represents NCBI reference of YP_004928144 with codon optimized for human cell expression

Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGAACTTTTCCGGGGAGTTCTCCAAGTTAGTTCCAATGTATTGGACTGCGCCAACGACAACACTGGTGGT
GTAGCCTTTTGGACCTGGACACATCAGATTGGGAGCCCTGACACATACAAATCGTTGATGGCTATTTA
CCTGTCCTCAGTAGCAAGCAAGCTGGACTTGACTGGCGGCCCTCTGGCAGGCTGCCTGTACTTCTTTTCAG
GCCGAGTGAATAAGTTTGAAGAGGGTTATCACATTCACGTGGTCATTGGTGGCCTGGTCTTAACCCGC
GCAATCTTACGGTGTGTGTTGAAGGGCTCTTCAATAACGTATTGTATCACTTCGTTACTGAAAATGTCAA
ATTGAAATTCCTGCCTGGCATGACCACAAAGGAAATACTTTCCGGATGGCGAGCAGTTTATCGAGAAC
TACCTGATGAAAAAGATTCCACTGAACGTGGTATGGTGCCTGACAAATATCGACGGCTATATCGACACAT
GTATATCAGCCACCTTCCGCAGAGGGGCTTGTACGCTAAGAAGCCCGGATCACCACCGCCATAAACGA
CACTAGCTCCGATGCCGCGAATCCTCCGGTACCGGAGCCGAGGTCGTGCCGTTTAAATGGGAAAGGCACT
AAAGCCAGCATCAAATTTAGACTATGGTTAATTGGCTTTGCGAGAACAGAGTTTTTACTGAAGATAAAT
GGAAGCTGGTCGACTTCAACCAGTACACATTGCTCAGCTCATCACATTCAGGATCCTTCAAATACAATC
AGCGCTGAAGCTGGCAATTTACAAGCGGACCAATCTCGTGCCGACAAGCACTTTCCTGTTGCACACAGAC
TTCGAGCAGGTAATGTGTATAAAGGATAACAAGATCGTCAAATGCTTCTGTGCCAGAATATGATCCCT
TGCTGGTGGGCCAACATGTGCTGAAGTGGATAGACAAGAAGTGGGTAACAAAAACACATTGTGGTTCTA
CGCCCCCGAGCACTGGTAAGACTAACCTCGCAATGGCAATCGCAAAAGTGTGCCTGTTTATGGAATG
GTCAATTTGGAACAACGAGAATTTCCATTCAATGACGTGGCCGGCAAGTCCCTTGTGGTGTGGGACGAGG
GGATCATCAAGAGTACCATTGTGGAAGCCGAAAAGCCATACTGGCGGCCAGCCTACGCGCGTGGATCA
GAAAATGCGCGGAGCGTCCGCCGTTCCCGGAGTCCCGTGTGAATTACAAGTAATGGAGATATTACTTTC
GTTGTGAGCGGGAATACCACAACAACGGTTCACGCCAAGGCGCTTAAGGAGCGAATGGTGAACCTTAACT
TCACAGTGCAGTGTTCACCGACATGGGGCTTCTGACAGAGGCAGACGTTTACGCAATGGTTGACATGGTG
TAATGCTCAGTCTGGGATCACTATGAGAATTGGGCTATCAACTACACCTTTGATTTCCAGGGATTAAT
GCTGATGCCCTTCATCCGGACCTCCAGACTACGCCTATTGTGACAGACACCAGCATCAGCAGCAGTGGCG
GGGAAAGTAGCGAAGAATTTCTGAGAGCTCTTTTTTCAACCTCATAACCCCGGGGGCATGGAACACGGA
GACTCCCAGGTCTAGCACTCCTATCCCAGGTACGTCTTCCGGCGAAAGCTTTGTCCGCTCCCGGTCAGC
AGCGAGGTGGTCGCCGTTTCATGGGAAGAAGCCTTCTACACGCCTCTGGCCGATCAGTTCAGGGAATGC
TGGTCCGTTGATGATTACGTCTGGGATGGGGTGGAGGGCCTGCCTGTCTGTTGTGTGCAGCACATCAATA
CTCCGGCGGAGGGCTGGGTCTGTGCCACACTGTATCAATGTAGGAGCGTGGTATAACGGATGGAAGTTT
AGGGAATTTACCCCGGACCTGGTGCATGCAGTGCACGTCGGGGCGTCCAATCCGTTTGTGATTGA
CCTGTAAGAAGTGCCTTACTTGTGAGGCTGCAATCCTTTGTGGACTATGAG

ACGCGTACGCGGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.</p>
Components:	<p>The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).</p>
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NC_000883.2 , YP_004928144
RefSeq ORF:	2013 bp
Locus ID:	1489597
MW:	74.1 kDa