

Product datasheet for **VC102101**

F (NC_001781) Virus Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	F (NC_001781) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	F Protein
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>The Viral ORF clone VC102101 represents NCBI reference of NP_056863 with codon optimized for human cell expression

Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGCTGCTTATTCACAGGCTTAGCGCTATTTTTTTGACCTTGCCATTAACGCCTTGTATCTGACAT
CCTCACAGAATATTACAGAAGAATTTTACCAGTCAACCTGCTCCGCGTTAGTAGAGGGTATTTCTCTGC
TCTCCGGACTGGCTGGTATACATCAGTCATCACGATTGAGTTGTCTAATATCAAGGAGACAAAATGCAAT
GGCACAGACACCAAAGTAACTCATCAAGCAGGAGTTGGACAAGTACAAAAATGCTGTGACGGAATTC
AACTGCTTATGCAGAATACCCCGCGCAACAATCGAGCAAGGAGGGAGGCTCCTCAGTACATGAATTA
TACAATTAATAACAACGAACCTTAATGTTTCAATCTCAAAAAAAGAAAAAGGAGATTCTCGGGTTC
CTGTTGGCGTGGCTCAGCCATTGCGAGCGGGATTGCTGTGTTCAAGTTCTCCACCTGGAGGGAGAAG
TGATAAGATTAAGAACGCGCTCCTGAGTACTAACAAGCCGTCGTTTCTCAGCAATGGGGTGTGAGT
GCTGACGAGTAAAGTTCTGGACCTGAAAACTATATCAACAATCAGCTGCTGCCATCGTCAACCAGCAG
AGCTGCCGATTAGCAACATTGAGACCGTAATCGAGTTCCAACAGAAGAACAGCCGGCTCCTGGAGATCA
ATCGAGAGTTCTCCGTTAATGCTGGTGTACTACGCCCTTAGCACTTACATGTTGACAAAATAGCGAGCT
CCTCTCCCTGATAAACGACATGCCGATTACTAACGATCAGAAAAAGCTGATGTCTAGTAACGTCCAGATC
GTCAGACAGCAGTCTTATTCAATAATGAGCATAATTAAGGAAGAAGTATTGGCCTACGTGGTCCAATGC
CAATCTACGGCGTGATCGATACCCCGTCTGGAACCTTCAACATAGCCCTCTGTGTACGACAAAATCAA
AGAGGGATCCAACATCTGTCTTACAGAACCGACCGGGCTGGTATTGCGACAACGCCGTAGCGCTCTCT
TTTTTCCACAGCCGATACCTGTAAGGTGCAGAGCAATAGGGTTTTTTGCGATACAATGAATCCCTGA
CGCTGCCCTCCGAAGTAAGCCTGTGCAATACGGATATCTCAATTCCAAATATGACTGTAATAAATGAC
ATCAAAAAACAGACATTTCTCCAGCGTGATCACCAGCCTCGGCGCAATTGTTAGTTGTTACGGCAAAACG
AAATGCACGGCTCCAATAAGAATAGAGGGATCATAAAACCTTCTCAAACGGCTGTGATTACGTAAGCA
ATAAGGGCGTGGACACCGTGAGCGTTGGAAATACTCTACTACGTCAATAAACTGGAGGGTAAAAACCT
CTATGTCAAGGGCGAGCCTATCATCAATTAATGATCCATTGGTGTTCCTAGTGACGAGTTCGATGCA
TCCATTTACAGGTTAACGAGAAGATTAATCAGTCCCTGGCATTATCCGGCGGAGCGACGAGCTGCTGC
ACAATGTTAATACCGTAAGTCCACTAACATCATGATAACCACAATTATCATTGTTATCATCGTCGT
CTTGCTGAGTCTATTGCTATCGGCCTTCTCTTTATTGCAAGGCCAAGAATACCCGGTACTCTCTCT
AAGGACCAACTGAGTGAATCAATAACATTGCATTCTCTAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>VC102101 representing NP_056863

Red=Cloning sites Green=Tags

MELLIHRLSAIFLTLAINALYLTSQQNITEEFYQSTCSAVSRGYFSALRTGWYTSVITIELSNIKETKCN
GTDTKVKLIKQELDKYKNAVTEQLLMQNTPAANNRARRREAPQYMYNTINTTKNLNVSISKKRKRRFLGF
LLGVGSAIASGIAVSKVLHLEGEVNIKIKNALLSTNKAVVSLSNGVSVLTSKVLDLKNYINNQLLPVYNQQ
SCRISNIETVIEFQQKNSRLLIENREFSVNAGVTTPLSTYMLTNSSELLSLINDMPITNDQKLMSSNVQI
VRQQSYSIMSIKEEVLAYVVQLPIYGVIDTPCWKLHTSPLCTTNIKEGSNICLTRDRGWYCDNAGSVS
FFPQADTKVQSNRVFCDTMNSLTLPSEVSLCNTDIFNSKYDCKIMTSKTDISSSVITSLGAIVSCYGKT
KCTASNKNRGIKTF SNGCDYVSNKGVDTVSVGNTLYYVKNLEGKNLYVKGEPIINYYDPLVFPDEFDA
SISQVNEKINQSLAFIRRSDELLHNVNTGKSTTNIMITTTIIVIIIVVLLSLIAIGLLLYKAKNTPVTL
KDQLSGINNIASF

TRTRRLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms:

https://cdn.origene.com/chromatograms/ja2124_e01.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NC_001781

ORF Size: 1722 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NC_001781.1](#), [NP_056863](#)

RefSeq ORF: 1722 bp
Locus ID: 1489825
MW: 63.7 kDa