

Product datasheet for **VC102169**

capsid protein VP3 (NC_010810) Virus Tagged ORF Clone

Product data:

| | |
|---------------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | capsid protein VP3 (NC_010810) Virus Tagged ORF Clone |
| Tag: | Myc-DDK |
| Symbol: | capsid protein VP3 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| ORF Nucleotide Sequence: | >The Viral ORF clone VC102169 represents NCBI reference of YP_001950224 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s) |

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCTCCAATCCCAGTTACCATTAGGGAACATGCTGGTACCTTCTATAGTACAACACCAGATACAACCTG
TGCCCATTTACGGAAAACTATCGCTACGCCTCCGACTATATGTGTGGAGAGTTTGATGACCTCCTGTC
CCTGTGTAAAATACCCACCTTCTGGCAACCCAACGACGGATCAAAGACTGCCATATATCTCAGCCAAC
AACACAGTCGGAGACCAGCCTCTGGCAACTTACCAGGTAACCTGAGTTGCAGCTGTATGGCTAATAGTA
TGCTCGCCCGCTCTCACGGAATTTAATCAATACAGAGGATCTCTGAATTATCTGTTCTGTTCCACCGG
CAGTGCTATGGCCAAGGGCAAATTTCTGATCAGTTACACTCCGCCTGGCGCTGGGAAACCTAAGGACAGG
GAACAGGCTATGCAAGCCACGTACGCGATTTGGGATTTGGGACTGAACAGCTCTTACAACCTTACTGTGC
CGTTTATCAGTCCCACACTACCGGCAGACCTCATAACCCACCCATCCATAACCAGTGTAGACGGTTG
GCTTAGTATTTGGCAACTCACTCCACTGACTTACCCTGCAAACACCCCAAGCGCCGACATCTTGACT
CTTGTCTCTGGCGGAGATGACTTCACGCTGCGAATGCCCGTTACACCATCCAAGTATACTCCTCAT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

RefSeq: NC_010810.1, YP_001950224
RefSeq ORF: 696 bp
MW: 25.5 kDa