

## Product datasheet for **VC102172**

### protein 2B (NC\_010810) Virus Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** protein 2B (NC\_010810) Virus Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** protein 2B  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >The Viral ORF clone VC102172 represents NCBI reference of YP\_001950227 with codon optimized for human cell expression  
**Red**=Cloning site **Blue**=ORF **Green**=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGGTGGCGTGCTCACTAAATCCCAGGCACCTATGTCTGGCCTCCAATCTATGCTGCTGCGGGCCATCG  
GCATTGAGGCAGACTGTACCGAATTCACCCGCGCTGTAACCTCATCACGGACTTGTGTAACACCTGGGA  
ATCCGCAAAGACAACCCCTCAGTTCTCTGAGTTCTGGACGAAGATGGTGATGCGGATAGTCAAGATGGTA  
GCTGCTTCGTTCTCTATCTCCACAACCCGATCTTACAACCACCTGTTTGTCTGCTCACTTATGGCTGGCA  
TCGACGTGCTGACCAACGACTCTGTTTTCAATTGGCTGAGCACCAAACCTGAGCAAGTTCTTTCATACCCC  
CGCACCTCCAATTGTGCCTTTGCTTCAACAGCAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >VC102172 representing YP\_001950227  
**Red**=Cloning sites **Green**=Tags

MGGVLTQSAPMSGLQSMLLRAIGIEADCTEFTRAVNLIITDLCNTWESAKTTLSPEFWTKMVMRIVKMW  
AASVLYLHNPDLTTTVCLSLMAGIDVLTNDSVFNWLSTKLSKFFHTPAPPVPLLQQQ

**TR**TRRLEQ**KL**I**SEEDLA**ND**ILDYK**DDDDK**V**

**Restriction Sites:** Sgfl-MluI



[View online »](#)

Cloning Scheme:



ACCN: NC\_010810

ORF Size: 384 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NC\\_010810.1](#), [YP\\_001950227](#)

RefSeq ORF: 384 bp

MW: 14.1 kDa