

## Product datasheet for **VC102183**

### 1C VP3 mature peptide (NC\_001489) Virus Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	1C VP3 mature peptide (NC_001489) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	1C VP3 mature peptide
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC102183 represents NCBI reference of NP_740551 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGATGATGCGGAATGAATTTAGGGTGAGCACTACGGAGAACGTGGTAAACCTGTCCAACATGAGGATG  
CACGGGCAAAGATGTCTTCGCCTTGGATCAGGAAGACTGGAAATCTGATCCCAGTCAGGGAGGGGGTAT  
CAAGATCACACACTTCACCACCTGGACTAGCATCCCCACCCTCGCCGCGCAGTTCCCATTCATGCTTCC  
GACTCCGTTGGACAGCAGATAAAAGTAATACCCGTGGACCCATATTTCTTCCAAATGACCAATACAAC  
CTGACCAAAAATGTATCACCGCCCTCGCCTCAATCTGCCAGATGTTTTGCTTTTGGCGGGGAGATCTCGT  
GTTGACTTCCAGGTGTTCCCCACGAAATATCATTCTGGCAGGCTGCTGTTTTGTTTCGTACCTGGTAAC  
GAACTTATCGACGTTAGCGGTATCACCCCTGAAACAAGCTACGACTGCCCTTGCGCCGTGATGGACATCA  
CCGGCGTACAGAGTACTCTGCGATTCCGAGTGCCTTGGATCAGTGACACACCGTACCGGGTAAATCGGTA  
CACCAATCCGCTCACCAGAAAGGTGAGTACACTGCTATTGGCAAGCTCATTGTATATTGTTACAACCGC  
TTGACTAGCCCTTCTAACGTCGCAAGCCATGTTGCGGTTAATGTGATTTGTCTGCCATAAACCTCGAGT  
GCTTCGCGCCCTGTACCATGCCATGGACGTGACAACCCAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >VC102183 representing NP\_740551  
 Red=Cloning sites Green=Tags

MMMRNEFRVSTTENVVNL SNYEDARAKMSFALDQEDWKS DPSQGGGK ITHFTTWT SIPTLAAQFPFNAS  
 DSVGQQIKVIPVDPYFFQMTNTNPDQKCITALASICQMF CFWRGDLVDFDQVFP TKYHSGRLLFCFVPGN  
 ELIDVSGITLKQATTAPCAVMDITGVQSTLRF RVPWISDTPYRVNRYTKSAHQKGEYTAIGKLIVYCYNR  
 LTSPSNVASHVRVNVYLSAINLECFAPLYHAMDVT TQ

TRTRRLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NC\_001489

**ORF Size:** 741 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NC\\_001489.1](#), [NP\\_740551](#)

**RefSeq ORF:** 741 bp

**MW:** 27.9 kDa