

## Product datasheet for **VC102190**

### 3C mature peptide (NC\_001489) Virus Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	3C mature peptide (NC_001489) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	3C mature peptide
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC102190 represents NCBI reference of NP_740558 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGTCTACTCTCGAGATTGCAGGACTCGTTCGAAAGAATCTGGTACAGTTCGGCGTCGGAGAGAAAAATG  
GATGTGTGAGATGGGTCATGAACGCTCTCGGGGTGAAAGATGACTGGCTGCTGGTTCCTTCTCATGCATA  
TAAGTTTGAGAAAGACTACGAAATGATGGAATTCTACTTTAACCGGGGGGCACATACTACAGTATTAGT  
GCTGGTAACGTAGTTATTCAGTCATTGATGTTGGATTTCAAGACGTGGTCTGATGAAAGTACCCACAA  
TTCCAAAAGTTCGGGATATAACTCAGCATTTTATCAAAAAGGGTATGTCGCCCGGGCTCTCAATAGACT  
GGCAACTCTCGTAACCACCGTGAATGGAACGCCCATGCTGATCTCTGAAGCCCCCTGAAGATGGAAGAG  
AAAGCCACATATGTGCACAAGAAGAATGACGGGACGACCGTGGACCTGACTGTCGACCAGGCATGGCGCG  
GCAAGGGCGAAGGCTCCTGGTATGTGTGGTGGCGCTCTGGTATCTAGTAACCAAAGTATCCAAAATGC  
CATCCTGGGCATCCAGTGGCCGGGGAAACAGTATACTCGTTGCGAACTGGTTACTCAGGAGATGTTCC  
CAGAATATTGATAAAAAAATCGAGTCTCAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >VC102190 representing NP\_740558  
 Red=Cloning sites Green=Tags

MSTLEIAGLVRKNLVQFGVGEKNGCVRWVMNALGVKDDWLLVPSHAYKFEKDYEMMEFYFNRRGGTYYSIS  
 AGNVVIQSLDVGFDVVLKMKVPTIPKFRDITQHFIKKGDVPRALNRLATLVTTVNGTPMLISEGPLKMEE  
 KATYVHKKNDGTTVDLTVDAQWRGKGEGLPGMCGGALVSSNQSISQNAILGIHVAGNSILVAKLVTQEMF  
 QNIDKKIESQ

TRTRRLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NC\_001489

**ORF Size:** 660 bp

**OTI Disclaimer:** The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NC\_001489.1, NP\_740558

RefSeq ORF: 660 bp

MW: 24.2 kDa