

## Product datasheet for VC102214

### protein F (NC\_004102) Virus Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	protein F (NC_004102) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	protein F
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC102214 represents NCBI reference of NP_803170 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGAGTACTAATCCGAAGCCACAGCGCAAACCTAACGTGACGCCAACCGTGGCGCATAGAACCAGCAGCT  
CTAGGGTGGCTGTTAGGTCACTTGTAGAATTCACCTGTTGTCGCGGGCGCACTTGACTGGGTGTGCGC  
ACGGAGGGGAAGGCTGCCTTCAGGAAGAAATCTGAGGTGGATGTCAGCCTCAGCCCTAGACACGTTGGC  
CCACGCGCAGGTCCAGGACTGAGTCCAGGCACCCTCGGACCCAGCATGGCGATGAGGGTTGCTGGCGGGC  
GAGATGGGTCTTGTCTCCCGTAGCCTTGGGACTCGCGGGTCTCCCAAACACCCGGTGTGGGCGCGC  
CATATGGGTACGATCTTCAATTCAGTGCAGCGCCGCTAGTCTACCTCTTGGGTACATACCGCAGCAGT  
GCACCACTGTTGGAAGCACTGCCCGCCGTGGCGAATGGCCTCTGGCTTTGGAAAACCGCC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:	>VC102214 representing NP_803170 Red=Cloning sites Green=Tags
-------------------	--

MSTNPKPQRKPNVTPTVAHRTSSSRVAVRSLVEFTCCRAGALDWVCARRGRLPSGRNLEVDVLSPRHVG  
PRAGPGLSPGTLGPMAMRVAGGRDGSCLPVALGLAGAPQTPGVGRAIWRSSIPLRAASPTSWGTYRSS  
APLLEALPGPWRMASGFWKTA

**TR**TRRLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms:	<a href="https://cdn.origene.com/chromatograms/mk8100_b10.zip">https://cdn.origene.com/chromatograms/mk8100_b10.zip</a>
----------------	---

Restriction Sites:	Sgfl-Mlul
--------------------	-----------



View online »

Cloning Scheme:



ACCN: NC\_004102

ORF Size: 483 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NC\\_004102.1, NP\\_803170](#)

RefSeq ORF: 483 bp

Locus ID: 951476

MW: 16.9 kDa