

Product datasheet for **VC102282**

HA (NC_006499) Virus Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	HA (NC_006499) Virus Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	HA
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>The Viral ORF clone VC102282 represents NCBI reference of YP_145800 with codon optimized for human cell expression Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCCCGGTTTCATCATCCTGTTCTGCTCCTGGCCCCGTGTACAGCAGACTGTGTCTGCGGAATCACC
CCGATACTACCTGGATCGGCGACAGCCGGTCCGACCAGTCCCGGGTGAACCAACAGAGCCTCGACCTGGT
GACCAACTTCAAGGGCATTCTGCAGGCCAAGAACGGCAACGGCCTGATGAAGCAGATGAGCGGCCGTTT
CCAGCGACTGGTACCAACCCACCACCAAGTACCGGATCCTCTACATCGGCACCAACGATTGTACCGAGG
GGCCCAACGACGTCATCATCCCCACCTCCATGACCCCTGGACAACGTGGCACGGGACCTCTACCTGGGCGC
CTGTCCGGGAGACGTGCGGGTGACCCCCACCTTCGTGGGCGCCGCGAGCTGGGCCTCATTGGGCGGACC
GACGCCCTGACCGAGTTCAGCGTGAAGGTGCTGACCTTCAACAACCCACCATTGTGGTGGTGGGCTCA
ATGGCATGTCCGGCATCTACAAGGTGTGCATCGCCGCCAGCAGCGCAACGTCCGGCGCGTGAACCTCGT
GAATGGGTGCGGCTACTTCAGCGCCCCCTGCGGTTTCGATAAATTCAAGGGCCAGATCTACGTGAGCGAC
ACCTTTGAGGTGCGGGGACCAAAAACAAGTGCCTGATCCTGCGGTCCAGCAGCAACGCCCCCTCTGCA
CCCACATCAAACGGAACATCGAGCTGGACGAGTACGTGGACACCCCAACACCGGCGGGTCTATCCCTC
CGACGGGTTTCAGAGCCTGCACGGCTCCGCCAGCATCCGGACCTTCCTCACCGAGGCCCTGACCTGCCCC
GGAGTGGATTGGGACCGGATCGACGCCCGCAGCTGCGAGTACGATTCCTGCCCAAGCTGGTCAAGGAGT
TCGACCAGACCGCCTGGGCAACACCGACCCAGATCATGCGGGAGCTCGAGGCCCAGAAGGAGATGAT
CGGCAAGCTCGGCCGGAACATCACCGACGTGAACAACCGGGTGGACGCCATCCCCCCCAGCTGAGCAAC
ATCTTCATCAGCATGGGGTGGCCGGCTTCGGAATCGCCCTGTTTCTGGCCGGCTGGAAGGCCTGCGTGT
GGATCGCTGCATTATGTACAAAAGCCGGGGCCGGAACCCCTGCCAACCTCAGCGTCGCA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >VC102282 representing YP_145800
Red=Cloning sites Green=Tags

MARFIIILFLLLAPVYSRLCLRNHPDTTWIGDSRSDQSRVNQQSLDLVTNFKGILQAKNGNGLMKQMSGRF
 PSDWYQPTTKYRILYIGTNDCTEGPNDVI IPTSMTLDNVARDL YLGACRGDVVRVPTTFVGAELGLIGRT
 DALTEFSVKVLTFNNPTIVVVGLNGMSGIYKVCIAASSGNVGGVNLVNGCGYFSAPLRFDNFKGQIYVSD
 TFEVRGBTKNKCIVILRSSNAPLCTHIKRNIELDEYVDTPTNTGGVYPSDGFDSLHGSASIRTFLEALTCP
 GVDWDRIDAASCEYDSCPKLVEKFDQTGLGNTDQIMRELEAQKEMIGKLRNITDVNNRVDAIPPQLSN
 IFISMGVAGFGIALFLAGWKACVWIAAFMYKSRGRNPPANLSVA

TRTRRLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NC_006499

ORF Size: 1182 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NC_006499.1](#), [YP_145800](#)

RefSeq ORF: 1182 bp

Locus ID: 3170814

MW: 43.1 kDa