

Product datasheet for **VC102288**

P7 (NC_006497) Virus Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: P7 (NC_006497) Virus Tagged ORF Clone
Tag: Myc-DDK
Symbol: P7
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >The Viral ORF clone VC102288 represents NCBI reference of YP_145797 with codon optimized for human cell expression
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**

ATGCGGGAGGCCAACCCCAAGCCCCGGGAGCTGATCCGGCACGCCCTGAAGAAGAAGAAGAGACCCGAGG
 TGGTCTACGCCATGGGCGTGTCTGCTGACACTGGGCGGCAGAGCGGCCTGACCGTGGAGTCCCAGTGCC
 CGAAGGCAAGACCGTGAAGGTCAAGACCCCTCAACCAGCTCGTCAACGGCATGATCAGCCGGCCACCATG
 ACTCTGTACTGCGTGAAGGACCCCGGCGGCATGGTACCCGATGCGGGACCACATCCGGA
 ACTGGCTGAAAGAGGAGCGGCTGCCAGGACCGCGCGGCGAGGAGAAATGGGCCATGGTGTACGG
 CATGATCAGCCCCGATATGGCCGAGGAAAAACCATGCTGAAGGAGCTGAAGACCATGCTCCATTCCCGG
 ATGCAGATGTACGCCCTGGGCGCCTCCAGCAAGGCCCTGGAGAACCTGGAGAAGGCCATCGTGGCCGCGAG
 TGCACCGGCTCCCGCCTCCTGCAGCACCGAGAAGATGGTGTCTCGGATACCTCAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >VC102288 representing YP_145797
Red=Cloning sites **Green**=Tags

MREANPKPRELIRHALKKKRPEVVYAMGVLLTLGGESGLTVEFPVPEGKTVKVKTLNQLVNGMISRATM
 TLYCVMKDPSSGMATLMRDHIRNWLKEESGCQDADGGEEKWAMVYGMISPDMAEEKTMLKELKTMHSR
 MQMYALGASSKALENLEKAIVAAVHRLPASCSTEKMLVLLGYLK



TRTRRLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



Cloning Scheme:

Cloning sites used for ORF Shutting:



* The last codon before the Stop codon of the ORF

ACCN: NC_006497

ORF Size: 549 bp

OTI Disclaimer: The molecular sequence of this clone can be viewed by clicking the "ORF Nucleotide Sequence" link above. This sequence represents the NCBI reference after codon optimization for human cell expression, and retaining the same decoded protein sequence. The stop codon in the native sequence was removed to create the in-frame c-terminal fusion with a Myc-DDK tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NC_006497.1](#), [YP_145797](#)

RefSeq ORF: 549 bp

Locus ID: 5075848

MW: 20.3 kDa